

VWR[®] TRACEABLE[®] DIGITAL HUMIDITY/TEMPERATURE/DEW POINT METER

Convenient triple memories recall MIN/MAX readouts for RH, Dew point, and Temperature

Individually calibrated

Individually certified—traceable to NIST

NEW



RH/Temperature/dew-point triple display

Convenient triple memories recall minimum/maximum readouts for relative humidity, dew point, and temperature. Simultaneously displays all three current readings. Always-on feature allows for monitoring conditions all night, over weekends, and for extended time periods. Oversized display supplies updates once a second.

SPECIFICATIONS

- Audible and visual LED alarm may be set to trigger at three different humidity readings (55%, 60% or 65%)
- Solid state sensors deliver response times of less than 30 seconds
- °F/°C switchable
- Size: 2.625 x 3.875 x 0.625 inches
- Weight: 6 ounces
- Supplied: wall mount, batteries, Traceable[®] Certificate

Cat. No.	Range	Resolution	Accuracy
	Humidity: 1 to 99%	0.1%	±3% mid-range, ±5% elsewhere
89500-398	Ambient Temperature: 14 to 140°F (-10 to 60°C)	0.1°	±1°C
	Ambient Dew point: 14 to 140°F (-10 to 60°C)	0.1°	na

Additional VWR® Traceable® Humidity Products

VWR® Traceable® Digital Humidity/Temperature Meter

Continuously/simultaneously displays relative humidity and temperature with dual-recall memories

- Ideal for routine measurements, monitoring round-the-clock, quality-control needs, and critical experimental requirements
- Minimum and maximum memories can be displayed or cleared at any time with the touch of a key
- Fulfills CLIA requirements, °F/°C switchable
- Supplied: wall mount, bench, stand, batteries, Traceable® Certificate



Cat. No.	Range	Resolution	Accuracy	Features
35519-047	Humidity: 25 to 95%RH	1%	±2% RH mid-range to ±4% RH elsewhere	MIN/MAX memory
	Temperature: 32 to 122°F (0 to 50°C)	1°	±1°C	

VWR® Traceable® Dew Point/Wet-Bulb/Humidity/Thermometer Alarm

Settable high/low alarms for all four readings

- Unit shows current minimum and maximum readings for dew point, wet-bulb, humidity and temperature
- Alarm sounds for 1 minute and red LED flashes until readings return to non-alarm conditions
- Supplied: stand, batteries, wall-mount, Traceable® Certificate



Cat. No.	Range	Resolution	Accuracy
36934-164	Humidity: 1 to 99%RH	0.1%	±4% between 20 to 80% RH
	Temperature: -40 to 158°F (-40 to 70°C)	0.1° below 100°	±1°C

VWR® Traceable® Memory-Card Humidity/Temperature/Dew Point

It records readings from once a minute up to once every 720 minutes on a SD memory card

- Datalogging capabilities that can capture 5.9 million readings with 256 mb removable SD memory card
- Visual and audible alarms indicate when temperature, humidity or dew point rises above and falls below user-defined set points
- USB memory card reader plugs directly into Windows and Mac—no software to install to retrieve date/time of day and readings of sensors
- Supplied: 256 mb SD memory card, USB memory card reader, magnetic strips, Velcro®, wall mount, bench stand, batteries, adaptor, Traceable® Certificate



Cat. No.	Range	Resolution	Accuracy	Dew point
89369-158	Humidity: 0.0 to 90.0%RH	0.1%	±5% RH at 25°C between 20 to 85% RH; ±7% RH elsewhere	-20.0 to 50.0°C with 0.1° resolution
	Temperature: 0 to 50.0°C	0.1°	±1°C	

Traceable to NIST—To assure accuracy an individually serial-numbered Traceable Certificate is provided from an ISO/IEC 17025:2005 (1750.01) calibration laboratory accredited by A2LA. It indicates traceability to standards provided by NIST (National Institute of Standards and Technology).



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.