



FEEL THE SW DIFFERENCE®


# DESIGNED TO ENHANCE WORKER PERFORMANCE.

## PowerForm®

**PF-95TL**  
POWDER-FREE ECOTEK SUSTAINABLE  
NITRILE EXAMINATION GLOVES\*

**VERSATILE SUSTAINABLE 5 MIL SINGLE-USE  
GLOVE COMBINES MAXIMUM COMFORT AND  
CONTROL WITH HIGH-QUALITY CONSTRUCTION**

- EcoTek® Sustainable Technology\*
- Silicone-free for critical environments
- Exceeds 1.5 AQL laboratory standards







Optimal thickness for durability and tactile sensitivity

Fully textured surface for wet and dry grip

Tear-resistant beaded cuff for easy donning

**PF** **PF-95TL**

WEIGHT	LENGTH	THICKNESS	FEATURES	COLOR
 6.5g (large size)	 9.5"	 5.0 mil	EcoTek® Clean Washed Air Inspected	 TEAL

**ECOTEK**




SW®  
33278 Central Avenue, Unit 102, Union City, CA 94587, USA  
Tel: +1.510.429.8692 | Fax: +1.510.487.5347  
Trademarks and registered trademarks are the property of SW and its affiliates.  
©2024 SW. All rights reserved.



# TECHNICAL DATA



**PF-95TL**  
**POWDER-FREE ECOTEK SUSTAINABLE**  
**NITRILE EXAMINATION GLOVES\***



## PRODUCT DESCRIPTION

### SPECIFICATIONS

<b>PRODUCT</b>	PowerForm®		
<b>GRADE</b>	Powder-Free Exam		
<b>MATERIAL</b>	Nitrile		
<b>THICKNESS</b>	<b>FINGER</b>	5.8mil	0.15mm
	<b>PALM</b>	5.0mil	0.13mm
	<b>CUFF</b>	3.9mil	0.10mm
<b>COLOR</b>	Teal		
<b>SURFACE</b>	Fully Textured		
<b>LENGTH</b>	9.5"	242mm	
<b>WEIGHT</b>	6.5g		
<b>FIT TYPE</b>	Ambidextrous		
<b>CUFF TYPE</b>	Beaded		
<b>WATERTIGHT AQL</b>	1.5		
<b>ORIGIN COUNTRY</b>	China		

### PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
<b>TENSILE</b>		
<b>STRENGTH:</b>	21 MPA	20 MPA
<b>ELONGATION:</b>	560%	540%

### INTEGRATED TECHNOLOGIES

**ECOTEK Sustainable Technology:** Product has been tested per ASTM D5526\*

100% Air Inspected


Clean Washed

\* Specifications based on size Large glove

## STANDARDS & CERTIFICATIONS

EN ISO 21420:2020  
 EN 455-2000-2015  
 ASTM D6319-19  
 ASTM D6978-05 (2019)  
 ASTM F1671-13  
 ASTM F739-12  
 Reach EU 1907/2006

EU 2016/425  
 CAT III Complex Design



Please visit [swssglobal.com/symbols](http://swssglobal.com/symbols) for information on standards used on this document.

## ORDERING

VWR SKU NO.	MFG NO.	SIZE	PALM WIDTH	PACKAGING
76051-116	N200362	S 6½ – 7	3.5" 89mm	100 gloves per box 10 boxes per case
76051-118	N200363	M 7½ – 8	3.8" 98mm	100 gloves per box 10 boxes per case
76051-120	N200364	L 8½ – 9	4.3" 108mm	100 gloves per box 10 boxes per case
76051-122	N200365	XL 9½ – 10	4.6" 116mm	100 gloves per box 10 boxes per case
76051-124	N200366	2XL 10½ – 11	4.9" 124mm	100 gloves per box 10 boxes per case

## APPLICATIONS

Chemical Handling

Keep out of sunlight. Store in a cool, dry place.  
 Keep away from sources of ozone or ignition.



R\_N20036X1121



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.

0524 Lit. No. 190219W

\* EcoTek gloves tested to ASTM D5526-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastic Materials Under Accelerated Landfill Conditions. Results are 92.6% biodegradation in 60% solids landfill at 945 days. Future results cannot be predicted/extrapolated. If you need more information on the test reports, send email to [cs@swssglobal.com](mailto:cs@swssglobal.com).