



G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

A Geno Technology, Inc. (USA) brand name

OneQuant™ HOOK™ Sulfo-NHS-SS-Biotin

(Cat. # 786-699)




think proteins! think G-Biosciences www.GBiosciences.com

INTRODUCTION

Sulfo-NHS-SS-Biotin is a commonly used amine reactive biotin-labeling reagent. OneQuant™ HOOK™ Sulfo-NHS-SS-Biotin are single aliquots of Sulfo-NHS-SS-Biotin reagent in convenient, premeasured microtubes. Each microtube contains 1mg of the biotin reagent. Simply dissolve one vial with appropriate solvent and add to the reaction.

SPECIFICATIONS

- **Formula:** C₁₉H₂₇O₉N₄S₄Na
- **Molecular Weight:** 606.69
- **Spacer Arm (Å):** 24.3
- **Reactive Group:** sulfo-NHS-ester
- **Membrane Permeable:** NO
- **Water Soluble:** YES
- **Cleavable/ Reversible:** YES
- **Reaction pH:** 7-9
- **Cas #:** 202057-28-1



ITEM(S) SUPPLIED (Cat. # 786-699)


Description	Size
OneQuant™ HOOK™ Sulfo-NHS-SS-Biotin	1mg/vial, 1 strip of 8 vials

STORAGE CONDITIONS

It is shipped at ambient temperature. Upon receipt store product at 4°C protected from moisture. Stable for one year.

RELATED PRODUCTS

Download our Protein Labeling and Conjugation Handbook.



<http://info.gbiosciences.com/complete-protein-labeling-conjugation-handbook>

For other related products, visit our website at www.GBiosciences.com or contact us.

Last saved: 9/25/2012 CMH



www.GBiosciences.com