

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Black Ink Mixture

SECTION 1: Identification of the substance/mixture and of the supplier

Product name: Black Ink Mixture

Manufacturer/Supplier Trade name:

Manufacturer/Supplier Article number: KEMBI1000-AA

Recommended uses of the product and restrictions on use: Laboratory

Manufacturer Details:

AquaPhoenix Scientific
9 Barnhart Drive, Hanover, PA 17331
(717) 632-1291

Supplier Details:

AquaPhoenix Scientific Inc.
9 Barnhart Drive, Hanover PA 17331
(717) 632-1291

Emergency telephone number:

Emergency Telephone No.: 800-255-3924

SECTION 2: Hazards identification

Classification of the substance or mixture: Not classified for physical or health hazards under GHS.

Signal word: None

Hazard statements: None

Precautionary statements:

If medical advice is needed have product container or label at hand.
Keep out of reach of children.
Read label before use.
Do not eat, drink or smoke when using this product.

Other Non-GHS Classification: None

SECTION 3: Composition/information on ingredients

Ingredients:

Ingredients:		
CAS N/A	Deep Blue Coloring	33.3 %
CAS N/A	Cherry Red Color, 1 Gallon	33.3 %
CAS 7732-18-5	Deionized Water	21.4 %
CAS N/A	Mordant Yellow 8	12 %
Percentages are by weight		

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Black Ink Mixture

SECTION 4: First aid measures

Description of first aid measures

After inhalation:

Move exposed individual to fresh air. Loosen clothing as necessary and position individual in a comfortable position. Seek medical advice if discomfort or irritation persists.

After skin contact:

Wash affected area with soap and water. Rinse thoroughly. Seek medical attention if irritation, discomfort or vomiting persists.

After eye contact:

Protect unexposed eye. Rinse/flush exposed eye(s) gently using water for 15-20 minutes. Remove contact lens(es) if able to do so during rinsing. Seek medical attention if irritation persists or if concerned.

After swallowing:

Rinse mouth thoroughly. Do not induce vomiting. Have exposed individual drink sips of water. Seek medical attention if irritation, discomfort or vomiting persists.

Most important symptoms and effects, both acute and delayed:

Irritation. Nausea. Headache. Shortness of breath.

Indication of any immediate medical attention and special treatment needed:

If seeking medical attention, provide SDS document to physician.

SECTION 5: Firefighting measures

Extinguishing media

Suitable extinguishing agents:

If in laboratory setting, follow laboratory fire suppression procedures. Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition.

Unsuitable extinguishing agents: None

Special hazards arising from the substance or mixture:

Combustion products may include carbon oxides or other toxic vapors.

Advice for firefighters:

Protective equipment: None

Additional information (precautions):

Move product containers away from fire or keep cool with water spray as a protective measure, where feasible.

SECTION 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures:

Wear protective equipment. Use respiratory protective device against the effects of fumes/dust/aerosol. Keep unprotected persons away. Ensure adequate ventilation. Keep away from ignition sources. Protect from heat. Stop the spill, if possible. Contain spilled material by diking or using inert absorbent. Transfer to a disposal or recovery container.

Environmental precautions:

Prevent from reaching drains, sewer or waterway. Collect contaminated soil for characterization per Section 13.

Methods and material for containment and cleaning up:

If in a laboratory setting, follow Chemical Hygiene Plan procedures. Collect liquids using vacuum or by use of absorbents. Place into properly labeled containers for recovery or disposal. If necessary, use trained response staff/contractor.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Black Ink Mixture

Reference to other sections: None

SECTION 7: Handling and storage

Precautions for safe handling:

Prevent formation of aerosols. Follow good hygiene procedures when handling chemical materials. Do not eat, drink, smoke, or use personal products when handling chemical substances. If in a laboratory setting, follow Chemical Hygiene Plan. Use only in well ventilated areas. Avoid splashes or spray in enclosed areas.

Conditions for safe storage, including any incompatibilities:

Store in a cool location. Provide ventilation for containers. Avoid storage near extreme heat, ignition sources or open flame. Store away from foodstuffs. Store away from oxidizing agents. Store in cool, dry conditions in well sealed containers. Keep container tightly closed.

SECTION 8: Exposure controls/personal protection



Control parameters:

No applicable occupational exposure limits.

Appropriate engineering controls:

Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use/handling. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of vapor or mists below the applicable workplace exposure limits (Occupational Exposure Limits-OELs) indicated above.

Respiratory protection:

Not required under normal conditions of use. Use suitable respiratory protective device when high concentrations are present. Use suitable respiratory protective device when aerosol or mist is formed. For spills, respiratory protection may be advisable.

Protection of skin:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation being used/handled. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.

Eye protection:

Safety glasses with side shields or goggles.

General hygienic measures:

The usual precautionary measures are to be adhered to when handling chemicals. Keep away from food, beverages and feed sources. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Do not inhale gases/fumes/dust/mist/vapor/aerosols. Avoid contact with the eyes and skin.

SECTION 9: Physical and chemical properties

Appearance (physical state, color):	Black liquid	Explosion limit lower: Explosion limit upper:	Not determined Not determined
Odor:	Odorless	Vapor pressure at 20°C:	Not determined
Odor threshold:	Not determined	Vapor density:	Not determined
pH-value:	Not determined	Relative density:	Not determined
Melting/Freezing point:	Not determined	Solubilities:	None
Boiling point/Boiling range:	Not determined	Partition coefficient (n-octanol/water):	Not determined

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Black Ink Mixture

Flash point (closed cup):	Not determined	Auto/Self-ignition temperature:	Not determined
Evaporation rate:	Not determined	Decomposition temperature:	Not determined
Flammability (solid, gaseous):	Not determined	Viscosity:	a. Kinematic: Not determined b. Dynamic: Not determined
Density at 20°C:	Not determined		

SECTION 10: Stability and reactivity

Reactivity: None

Chemical stability:

No decomposition if used and stored according to specifications.

Possible hazardous reactions: None

Conditions to avoid:

Store away from oxidizing agents, strong acids or bases.

Incompatible materials:

Strong acids. Strong bases.

Hazardous decomposition products:

Carbon oxides (CO, CO₂).

SECTION 11: Toxicological information

Acute Toxicity: None

Chronic Toxicity: No additional information.

Skin corrosion/irritation: No additional information.

Serious eye damage/irritation: No additional information.

Respiratory or skin sensitization: No additional information.

Carcinogenicity: No additional information.

Germ cell mutagenicity: No additional information.

Reproductive Toxicity: No additional information.

STOT-single and repeated exposure: No additional information.

Additional toxicological information: No additional information.

SECTION 12: Ecological information

Ecotoxicity: No additional information.

Persistence and degradability:

Readily degradable in the environment.

Bioaccumulative potential: No additional information.

Mobility in soil:

Aqueous solution has high mobility in soil.

Other adverse effects: No additional information.

SECTION 13: Disposal considerations

Waste disposal recommendations:

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Black Ink Mixture

Product/containers must not be disposed together with household garbage. Do not allow product to reach sewage system or open water. It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities (US 40CFR262.11). Consult federal state/ provincial and local regulations regarding the proper disposal of waste material that may incorporate some amount of this product.

SECTION 14: Transport information

US DOT

UN Number:

ADR, ADN, DOT, IMDG, IATA

Not Regulated.

Limited Quantity Exception:

None

Bulk:

RQ (if applicable): None

Proper shipping Name: Not Regulated.

Hazard Class: None

Packing Group: Not Regulated.

Marine Pollutant (if applicable): No additional information.

Comments: None

Non Bulk:

RQ (if applicable): None

Proper shipping Name: Not Regulated.

Hazard Class: None

Packing Group: Not Regulated.

Marine Pollutant (if applicable): No additional information.

Comments: None

SECTION 15: Regulatory information

United States (USA)

SARA Section 311/312 (Specific toxic chemical listings):

None of the ingredients are listed.

SARA Section 313 (Specific toxic chemical listings):

None of the ingredients are listed.

RCRA (hazardous waste code):

None of the ingredients are listed.

TSCA (Toxic Substances Control Act) :

All ingredients are listed.

CERCLA (Comprehensive Environmental Response, Compensation, and Liability Act):

None of the ingredients are listed.

Proposition 65 (California):

Chemicals known to cause cancer:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients are listed.

Chemicals known to cause developmental toxicity:

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Black Ink Mixture

None of the ingredients are listed.

Canada

Canadian Domestic Substances List (DSL) :

All ingredients are listed.

SECTION 16: Other information

This product has been classified in accordance with hazard criteria of the Controlled Products Regulations and the SDS contains all the information required by the Controlled Products Regulations. Note. The responsibility to provide a safe workplace remains with the user. The user should consider the health hazards and safety information contained herein as a guide and should take those precautions required in an individual operation to instruct employees and develop work practice procedures for a safe work environment. The information contained herein is, to the best of our knowledge and belief, accurate. However, since the conditions of handling and use are beyond our control, we make no guarantee of results, and assume no liability for damages incurred by the use of this material. It is the responsibility of the user to comply with all applicable laws and regulations applicable to this material.

NFPA: 0-0-0

HMIS: 0-0-0

GHS Full Text Phrases: None

Abbreviations and Acronyms:

IMDG International Maritime Code for Dangerous Goods.
IATA International Air Transport Association.
GHS Globally Harmonized System of Classification and Labelling of Chemicals.
ACGIH American Conference of Governmental Industrial Hygienists.
CAS Chemical Abstracts Service (division of the American Chemical Society).
NFPA National Fire Protection Association (USA).
HMIS Hazardous Materials Identification System (USA).
WHMIS Workplace Hazardous Materials Information System (Canada).
DNEL Derived No-Effect Level (REACH).
PNEC. Predicted No-Effect Concentration (REACH).
CFR Code of Federal Regulations (USA).
SARA Superfund Amendments and Reauthorization Act (USA).
RCRA. Resource Conservation and Recovery Act (USA).
TSCA. Toxic Substances Control Act (USA).
NPRI National Pollutant Release Inventory (Canada).
DOT US Department of Transportation.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder

SECTION 1: Identification of the substance/mixture and of the supplier

Product name: Sulfur Powder

Manufacturer/Supplier Trade name:

Manufacturer/Supplier Article number: KEMSU1500-12G

Recommended uses of the product and restrictions on use:

Manufacturer Details:

AquaPhoenix Scientific, Inc.
9 Barnhart Drive
Hanover, PA 17331
1-717-632-1291

Supplier Details:

AquaPhoenix Scientific Inc.
9 Barnhart Drive, Hanover PA 17331
(717) 632-1291

Emergency telephone number:

ChemTel: (24-hour) (US and Canada)
1-(800)-255-3924

SECTION 2: Hazards identification

Classification of the substance or mixture:



Irritant

Skin irritation, category 2



Flammable

Flammable solids, category 2

Flammable Solid 2.

Skin corrosion/irritation - Category 2.

HNOC: Combustible Dust.

Signal word: Warning

Hazard statements:

Causes skin irritation.

Precautionary statements:

If medical advice is needed have product container or label at hand.

Keep out of reach of children.

Read label before use.

Do not eat, drink or smoke when using this product.

Wash skin thoroughly after handling.

Wear protective gloves/protective clothing/eye protection/face protection.

Keep away from heat/sparks/open flames/hot surfaces. No smoking.

Ground/bond container and receiving equipment.

Use explosion-proof electrical/ventilating/light/equipment.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder

IF ON SKIN: Wash with soap and water.
If skin irritation occurs: Get medical advice/attention.
Remove/Take off immediately all contaminated clothing.
In case of fire, use agents recommended in section 5 for extinction.
Dispose of contents and container as instructed in Section 13.

Other Non-GHS Classification: None

SECTION 3: Composition/information on ingredients

Ingredients:

Ingredients:		
CAS 7704-34-9	Sulfur	>99 %
Percentages are by weight		

SECTION 4: First aid measures

Description of first aid measures

After inhalation:

Loosen clothing as necessary and position individual in a comfortable position. Move exposed to fresh air. Give artificial respiration if necessary. If breathing is difficult give oxygen. Get medical assistance if cough or other symptoms appear.

After skin contact:

Rinse/flush exposed skin gently using soap and water for 15-20 minutes. Seek medical advice if discomfort or irritation persists.

After eye contact:

Protect unexposed eye. Rinse/flush exposed eye(s) gently using water for 15-20 minutes. Remove contact lens(es) if able to do so during rinsing. Seek medical attention if irritation persists or if concerned.

After swallowing:

Rinse mouth thoroughly. Do not induce vomiting. Have exposed individual drink sips of water. Seek medical attention if irritation, discomfort or vomiting persists. Never give anything by mouth to an unconscious person.

Most important symptoms and effects, both acute and delayed:

Can produce delayed pulmonary edema. May cause central nervous system effects. Nausea. Headache. Shortness of breath. Olfactory fatigue may occur. Irritation- all routes of exposure. Diarrhea. Vomiting. Prolonged or repeated skin contact may cause dermatitis. Chronic inhalation may cause effects similar to those of acute inhalation.

Indication of any immediate medical attention and special treatment needed:

If seeking medical attention, provide SDS document to physician. Physician should treat symptomatically.

SECTION 5: Firefighting measures

Extinguishing media

Suitable extinguishing agents:

Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition. Water spray may cause frothing. For large fires, use water spray, fog or regular foam. Contact professional fire-fighters immediately. For small fires, use dry chemical, carbon dioxide, sand, earth, water spray or regular foam. Cool containers with water.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder

Unsuitable extinguishing agents: None

Special hazards arising from the substance or mixture:

Thermal decomposition can lead to release of irritating gases and vapors. sulfur oxides. During a fire, irritating and highly toxic gases may be generated by thermal decomposition or combustion. Dust can be an explosion hazard when exposed to heat or flame. This material in sufficient quantity and reduced particle size is capable of creating a dust explosion. Flammable solid. May burn rapidly with flare burning effect. May re-ignite after fire is extinguished.

Advice for firefighters:

Protective equipment:

Use NIOSH-approved respiratory protection/breathing apparatus. Evacuate area and fight fire from a safe distance. As in any fire, wear self-contained breathing apparatus in pressure-demand, MSHA/NIOSH (approved or equivalent), and full protective gear.

Additional information (precautions):

Move product containers away from fire or keep cool with water spray as a protective measure, where feasible. Use spark-proof tools and explosion-proof equipment. Avoid generating dust; fine dust dispersed in air in sufficient concentrations, and in the presence of an ignition source is a potential dust explosion hazard. Avoid inhaling gases, fumes, dust, mist, vapor, and aerosols. Avoid contact with skin, eyes, and clothing.

SECTION 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures:

Wear protective equipment. Use spark-proof tools and explosion-proof equipment. Ensure that air-handling systems are operational. Ensure adequate ventilation.

Environmental precautions:

Prevent from reaching drains, sewer or waterway. Collect contaminated soil for characterization per Section 13. Should not be released into environment.

Methods and material for containment and cleaning up:

Keep in suitable closed containers for disposal. Wear protective eyewear, gloves, and clothing. Always obey local regulations. Vacuum or sweep up material and place into a suitable disposal container. Clean up spills immediately, observing precautions in Section 8. Scoop up with a non-sparking tool, then place into a suitable container for disposal. Avoid generating dust. Remove all sources of ignition. Provide ventilation. Refer to Section 8. Dust deposits should not be allowed to accumulate on surfaces, as these may form an explosive mixture if they are released into the atmosphere in sufficient concentration. Avoid dispersal of dust in the air. Collect solids in powder form using vacuum with HEPA filter. Evacuate personnel to safe areas.

Reference to other sections: None

SECTION 7: Handling and storage

Precautions for safe handling:

Follow good hygiene procedures when handling chemical materials. Refer to Section 8. Wash thoroughly after handling. Remove contaminated clothing and wash before reuse. Minimize dust generation and accumulation. May form flammable dust-air mixtures. Avoid contact with eyes, skin, and clothing. Empty containers retain product residue, (liquid and/or vapor), and can be dangerous. Keep container tightly closed. Use with adequate ventilation. Keep away from heat and sources of ignition. Do not eat, drink, smoke, or use personal products when handling chemical substances. Avoid contact with eyes, skin, and clothing.

Conditions for safe storage, including any incompatibilities:

Protect from freezing and physical damage. Keep away from food and beverages. Keep away from heat, sparks, and flame. Keep away from sources of ignition. Store in a tightly closed container. Store in a cool, dry, well-ventilated area away from incompatible substances. Store with like hazards.

SECTION 8: Exposure controls/personal protection

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder



Control parameters:	No applicable occupational exposure limits.
Appropriate engineering controls:	Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use/handling. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of vapor or dusts (total/respirable) below the applicable workplace exposure limits (Occupational Exposure Limits-OELs) indicated above. It is recommended that all dust control equipment such as local exhaust ventilation and material transport systems involved in handling of this product contain explosion relief vents or an explosion suppression system or an oxygen deficient environment. Ensure that dust-handling systems (such as exhaust ducts, dust collectors, vessels, and processing equipment) are designed in a manner to prevent the escape of dust into the work area (i.e., there is no leakage from the equipment). Use under a fume hood.
Respiratory protection:	Not required under normal conditions of use. Where risk assessment shows air-purifying respirators are appropriate use a full-face particle respirator type N100 (US) or type P3 (EN 143) respirator cartridges as a backup to engineering controls. When necessary use NIOSH approved breathing equipment.
Protection of skin:	Select glove material impermeable and resistant to the substance. Select glove material based on rates of diffusion and degradation. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Use proper glove removal technique without touching outer surface. Avoid skin contact with used gloves. Wear protective clothing.
Eye protection:	Wear equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU). Safety glasses or goggles are appropriate eye protection.
General hygienic measures:	Perform routine housekeeping. Wash hands before breaks and at the end of work. Avoid contact with skin, eyes, and clothing. Before wearing wash contaminated clothing.

SECTION 9: Physical and chemical properties

Appearance (physical state, color):	Yellow solid	Explosion limit lower: Explosion limit upper:	3.3 %(V) 46 %(V)
Odor:	Rotten egg - like	Vapor pressure at 20°C:	10 hPa (8 mmHg) at 246 °C (475 °F) 1 hPa (1 mmHg) at 183.8 °C (362.8 °F)
Odor threshold:	Not determined	Vapor density:	Not determined
pH-value:	Not determined	Relative density:	Not determined
Melting/Freezing point:	Melting point/range:	Solubilities:	Insoluble in water
Boiling point/Boiling range:	444.7 °C (832.5 °F)	Partition coefficient (n-octanol/water):	Not determined
Flash point (closed cup):	207 °C (405 °F)	Auto/Self-ignition temperature:	232 °C (450 °F)
Evaporation rate:	Not determined	Decomposition temperature:	Not determined

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder

Flammability (solid, gaseous):	May form combustible dust concentrations in air.	Viscosity:	a. Kinematic: 8 mm ² /s at 140 °C (284 °F) - b. Dynamic: Not determined
Density at 20°C:	Not determined		
Specific Gravity:	2.07		

SECTION 10: Stability and reactivity

Reactivity:

Nonreactive under normal conditions.

Chemical stability:

Stable under normal conditions.

Possible hazardous reactions:

None under normal processing.

Conditions to avoid:

Incompatible Materials.

Incompatible materials:

Strong acids. Strong bases. Oxidizing agents.

Hazardous decomposition products:

Sulfur oxides (SO_x), including sulfur oxide and sulfur dioxide.

SECTION 11: Toxicological information

Acute Toxicity:

Dermal:

LD50 Dermal - Rabbit - > 2,000 mg/kg.

Chronic Toxicity: No additional information.

Skin corrosion/irritation: No additional information.

Serious eye damage/irritation: No additional information.

Respiratory or skin sensitization: No additional information.

Carcinogenicity:

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

Germ cell mutagenicity: No additional information.

Reproductive Toxicity: No additional information.

STOT-single and repeated exposure: No additional information.

Additional toxicological information: No additional information.

SECTION 12: Ecological information

Ecotoxicity:

Toxicity to fish , LC50 - Oncorhynchus mykiss (rainbow trout) - > 180 mg/l - 96 h LC50 - other fish - 866 mg/l - 96 h.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder

Toxicity to daphnia and other aquatic invertebrates, EC50 - Daphnia magna (Water flea) - > 5,000 mg/l - 48 h.

Persistence and degradability: No additional information.

Bioaccumulative potential: No additional information.

Mobility in soil: No additional information.

Other adverse effects: No additional information.

SECTION 13: Disposal considerations

Waste disposal recommendations:

Contact a licensed professional waste disposal service to dispose of this material. Dispose of empty containers as unused product. Product or containers must not be disposed with household garbage. It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities (US 40CFR262.11). Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. Chemical waste generators must also consult local, regional, and national hazardous waste regulations. Ensure complete and accurate classification.

SECTION 14: Transport information

US DOT

UN Number:

ADR, ADN, DOT, IMDG, IATA 1350

Limited Quantity Exception: None

Bulk:

RQ (if applicable): None

Proper shipping Name: Sulfur.

Hazard Class: 4

Packing Group: III.

Marine Pollutant (if applicable): No additional information.

Comments: None

Non Bulk:

RQ (if applicable): None

Proper shipping Name: Sulfur.

Hazard Class: 4

Packing Group: III.

Marine Pollutant (if applicable): No additional information.

Comments: None



SECTION 15: Regulatory information

United States (USA)

SARA Section 311/312 (Specific toxic chemical listings):

Acute, Fire

SARA Section 313 (Specific toxic chemical listings):

None of the ingredients are listed.

RCRA (hazardous waste code):

None of the ingredients are listed.

TSCA (Toxic Substances Control Act) :

All ingredients are listed.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder

CERCLA (Comprehensive Environmental Response, Compensation, and Liability Act):

None of the ingredients are listed.

Proposition 65 (California):

Chemicals known to cause cancer:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients are listed.

Chemicals known to cause developmental toxicity:

None of the ingredients are listed.

Canada

Canadian Domestic Substances List (DSL) :

All ingredients are listed.

SECTION 16: Other information

This product has been classified in accordance with hazard criteria of the Controlled Products Regulations and the SDS contains all the information required by the Controlled Products Regulations. Note. The responsibility to provide a safe workplace remains with the user. The user should consider the health hazards and safety information contained herein as a guide and should take those precautions required in an individual operation to instruct employees and develop work practice procedures for a safe work environment. The information contained herein is, to the best of our knowledge and belief, accurate. However, since the conditions of handling and use are beyond our control, we make no guarantee of results, and assume no liability for damages incurred by the use of this material. It is the responsibility of the user to comply with all applicable laws and regulations applicable to this material.

NFPA: 2-0-2

HMIS: 2-0-2

GHS Full Text Phrases: None

Abbreviations and Acronyms:

IMDG International Maritime Code for Dangerous Goods.
IATA International Air Transport Association.
GHS Globally Harmonized System of Classification and Labelling of Chemicals.
ACGIH American Conference of Governmental Industrial Hygienists.
CAS Chemical Abstracts Service (division of the American Chemical Society).
NFPA National Fire Protection Association (USA).
HMIS Hazardous Materials Identification System (USA).
WHMIS Workplace Hazardous Materials Information System (Canada).
DNEL Derived No-Effect Level (REACH).
PNEC. Predicted No-Effect Concentration (REACH).
CFR Code of Federal Regulations (USA).
SARA Superfund Amendments and Reauthorization Act (USA).
RCRA. Resource Conservation and Recovery Act (USA).
TSCA. Toxic Substances Control Act (USA).
NPRI National Pollutant Release Inventory (Canada).

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 10.24.2014

Sulfur Powder

DOT US Department of Transportation.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.05.2014

Sodium Chloride, ACS Grade

SECTION 1: Identification of the substance/mixture and of the supplier

Product name: Sodium Chloride, ACS Grade

Manufacturer/Supplier Trade name:

Manufacturer/Supplier Article number: KEMSC1010-12G

Recommended uses of the product and restrictions on use:

Manufacturer Details:

AquaPhoenix Scientific, Inc.
9 Barnhart Drive
Hanover, PA 17331
1-717-632-1291

Supplier Details:

AquaPhoenix Scientific Inc.
9 Barnhart Drive, Hanover PA 17331
(717) 632-1291

Emergency telephone number:

ChemTel: (24-hour) (US and Canada)
1-(800)-255-3924

SECTION 2: Hazards identification

Classification of the substance or mixture:



Irritant

Eye irritation, category 2B

Eye Irritation 2.

Signal word: Warning

Hazard statements:

Causes serious eye irritation.

Precautionary statements:

If medical advice is needed have product container or label at hand.

Keep out of reach of children.

Read label before use.

Wear protective gloves/protective clothing/eye protection/face protection.

Do not eat, drink or smoke when using this product.

Wash skin thoroughly after handling.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do - continue rinsing.

If eye irritation persists get medical advice/attention.

Other Non-GHS Classification: None

SECTION 3: Composition/information on ingredients

Ingredients:

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.05.2014

Sodium Chloride, ACS Grade

Ingredients:

CAS 7647-14-5	Sodium Chloride, ACS	100 %
Percentages are by weight		

SECTION 4: First aid measures

Description of first aid measures

After inhalation:

Move exposed individual to fresh air. Loosen clothing as necessary and position individual in a comfortable position. Seek medical advice if discomfort or irritation persists. Remove to fresh air. Give artificial respiration if necessary. If breathing is difficult, give oxygen.

After skin contact:

Wash affected area with soap and water. Wash affected area with soap and water. Rinse thoroughly. Seek medical attention if irritation persists or if concerned. Flush with water for 15 minutes. Rinse thoroughly. Seek medical attention if symptoms develop or persist.

After eye contact:

Protect unexposed eye. Remove contact lens(es) if able to do so during rinsing. Seek medical attention if irritation persists or if concerned. Immediately flush eyes with water for at least 15 minutes. Immediately get medical assistance.

After swallowing:

Rinse mouth thoroughly. Do not induce vomiting. Seek medical assistance. Dilute with water or milk.

Most important symptoms and effects, both acute and delayed:

Nausea. Headache. Shortness of breath. Irritation- all routes of exposure.

Indication of any immediate medical attention and special treatment needed:

If seeking medical attention, provide SDS document to physician. Note to physician: Treat symptomatically.

SECTION 5: Firefighting measures

Extinguishing media

Suitable extinguishing agents:

If in laboratory setting, follow laboratory fire suppression procedures. Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition.

Unsuitable extinguishing agents: None

Special hazards arising from the substance or mixture:

Combustion products may include sodium oxides or other toxic vapors. Avoid generating dust; fine dust dispersed in air in sufficient concentrations, and in the presence of an ignition source is a potential dust explosion hazard.

Advice for firefighters:

Protective equipment:

Use NIOSH-approved breathing equipment.

Additional information (precautions):

Move product containers away from fire or keep cool with water spray as a protective measure, where feasible. Use non-sparking equipment/tools.

SECTION 6: Accidental release measures

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.05.2014

Sodium Chloride, ACS Grade

Personal precautions, protective equipment and emergency procedures:

Wear protective equipment. Avoid contact with skin and eyes. Always obey local regulations. Avoid contact skin, eyes, and clothing. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed. Use respiratory protective device against the effects of fumes/dust/aerosol. Keep unprotected persons away. Ensure adequate ventilation. Keep away from ignition sources. Protect from heat. Stop the spill, if possible. Contain spilled material by diking or using inert absorbent. Transfer to a disposal or recovery container.

Environmental precautions:

Prevent from reaching drains, sewer or waterway. Collect contaminated soil for characterization per Section 13. Should not be released into the environment.

Methods and material for containment and cleaning up:

If in a laboratory setting, follow Chemical Hygiene Plan procedures. Collect liquids using vacuum or by use of absorbents. Place into properly labeled containers for recovery or disposal. If necessary, use trained response staff/contractor. Dust deposits should not be allowed to accumulate on surfaces, as these may form an explosive mixture if they are released into the atmosphere in sufficient concentration. Avoid dispersal of dust in the air. Always obey local regulations.

Reference to other sections: None

SECTION 7: Handling and storage

Precautions for safe handling:

Follow good hygiene procedures when handling chemical materials. Protect from freezing and physical damage. Keep away from sources of ignition. Store protected from moisture and direct sunlight. Wash hands after handling. Do not eat, drink, smoke, or use personal products when handling chemical substances. If in a laboratory setting, follow Chemical Hygiene Plan. Use only in well ventilated areas. Avoid generation of dust or fine particulate. Avoid dispersal of dust in the air. Routine housekeeping should be instituted to ensure that dusts do not accumulate on surfaces. Dry powders can build static electricity charges when subjected to the friction of transfer and mixing operations. Avoid contact with eyes, skin, and clothing.

Conditions for safe storage, including any incompatibilities:

Store in a cool location. Provide ventilation for containers. Avoid storage near extreme heat, ignition sources or open flame. Store away from foodstuffs. Store away from oxidizing agents. Store in cool, dry conditions in well sealed containers. Keep container tightly closed. Protect from freezing and physical damage. Keep away from sources of ignition. Store protected from moisture and direct sunlight.

SECTION 8: Exposure controls/personal protection



Control parameters:

7647-14-5, ACGIH TLV TWA (inhalable particles), 10 mg/m³.
7647-14-5, OSHA PEL TWA (Total Dust), 15 mg/m³ (50 mppcf*).
, , *mppcf = Millions of particles per cubic foot of air.

Appropriate engineering controls:

Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use/handling. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of vapor or dusts (total/respirable) below the applicable workplace exposure limits (Occupational Exposure Limits-OELs) indicated above. Ensure that dust-handling systems (such as exhaust ducts, dust collectors, vessels, and processing equipment) are designed in a manner to prevent the escape of dust into the work area (i.e., there is no leakage from the equipment).

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.05.2014

Sodium Chloride, ACS Grade

Respiratory protection:	Not required under normal conditions of use. Use suitable respiratory protective device when high concentrations are present. Use suitable respiratory protective device when aerosol or mist is formed. For spills, respiratory protection may be advisable.
Protection of skin:	The glove material has to be impermeable and resistant to the product/ the substance/ the preparation being used/handled. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation. Wear protective clothing.
Eye protection:	Safety glasses with side shields or goggles.
General hygienic measures:	The usual precautionary measures are to be adhered to when handling chemicals. Keep away from food, beverages and feed sources. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Do not inhale gases/fumes/dust/mist/vapor/aerosols. Avoid contact with the eyes and skin.

SECTION 9: Physical and chemical properties

Appearance (physical state, color):	White solid	Explosion limit lower: Explosion limit upper:	Not determined Not determined
Odor:	Odorless	Vapor pressure at 20°C:	1 mmHg @ 865°C
Odor threshold:	Not determined	Vapor density:	>1
pH-value:	Not determined	Relative density:	Not determined
Melting/Freezing point:	801°C	Solubilities:	Soluble in water.
Boiling point/Boiling range:	1461°C	Partition coefficient (n-octanol/water):	Not determined
Flash point (closed cup):	Not determined	Auto/Self-ignition temperature:	Not determined
Evaporation rate:	Not determined	Decomposition temperature:	Not determined
Flammability (solid, gaseous):	Not determined	Viscosity:	a. Kinematic: Not determined b. Dynamic: Not determined
Density at 20°C:	Not determined		
Specific Gravity	2.165		
Molecular Weight:	58.44 g/mol		

SECTION 10: Stability and reactivity

Reactivity:

Material is hygroscopic.

Chemical stability:

No decomposition if used and stored according to specifications.

Possible hazardous reactions:

None under normal processing.

Conditions to avoid:

Store away from oxidizing agents, strong acids or bases. Incompatible materials, dust generation, combustible materials, exposure to moist air or water. Excess heat.

Incompatible materials:

Metals. Strong oxidizers. Strong acids. Strong bases. Incompatible materials, dust generation, combustible

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.05.2014

Sodium Chloride, ACS Grade

materials, exposure to moist air or water.

Hazardous decomposition products:

May evolve chlorine gas when in contact with strong acids. Sodium/sodium oxides. Hydrogen chloride gas.

SECTION 11: Toxicological information

Acute Toxicity:

Dermal:

LD50 dermal-rabbit (7647-14-5) > 10gm/kg.

Chronic Toxicity: No additional information.

Skin corrosion/irritation: No additional information.

Serious eye damage/irritation: No additional information.

Respiratory or skin sensitization: No additional information.

Carcinogenicity: No additional information.

Germ cell mutagenicity: No additional information.

Reproductive Toxicity: No additional information.

STOT-single and repeated exposure: No additional information.

Additional toxicological information: No additional information.

SECTION 12: Ecological information

Ecotoxicity:

Fish (acute 7647-14-5), 96 Hr LC50 *Lepomis macrochirus*: 5560 - 6080 mg/L [flow-through]; 96 Hr LC50 *Lepomis macrochirus*: 12946 mg/L [static]; 96 Hr LC50 *Pimephales promelas*: 6020 - 7070 mg/L [static]; 96 Hr LC50 *Pimephales promelas*: 7050 mg/L [semi-static]; 96 Hr LC50 *Pimephales promelas*: 6420 - 6700 mg/L [static]; 96 Hr LC50 *Oncorhynchus mykiss*: 4747 - 7824 mg/L [flow-through].

Crustacea (acute 7647-14-5): , 48 Hr EC50 *Daphnia magna*: 1000 mg/L; 48 Hr EC50 *Daphnia magna*: 340.7 - 469.2 mg/L [Static].

Persistence and degradability:

Can attenuate over time. Large amounts can persist in the environment.

Bioaccumulative potential: No additional information.

Mobility in soil:

Soluble in water; thus mobile along soil/water interface.

Other adverse effects:

Should not be released into environment.

SECTION 13: Disposal considerations

Waste disposal recommendations:

Product/containers must not be disposed together with household garbage. Do not allow product to reach sewage system or open water. It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities (US 40CFR262.11). Consult federal state/ provincial and local regulations regarding the proper disposal of waste material that may incorporate some amount of this product.

SECTION 14: Transport information

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.05.2014

Sodium Chloride, ACS Grade

US DOT

UN Number:

ADR, ADN, DOT, IMDG, IATA

Not Regulated.

Limited Quantity Exception:

None

Bulk:

RQ (if applicable): None

Proper shipping Name: Not Regulated.

Hazard Class: None

Packing Group: Not Regulated.

Marine Pollutant (if applicable): No additional information.

Comments: None

Non Bulk:

RQ (if applicable): None

Proper shipping Name: Not Regulated.

Hazard Class: None

Packing Group: Not Regulated.

Marine Pollutant (if applicable): No additional information.

Comments: None

SECTION 15: Regulatory information

United States (USA)

SARA Section 311/312 (Specific toxic chemical listings):

Acute

SARA Section 313 (Specific toxic chemical listings):

None of the ingredients are listed.

RCRA (hazardous waste code):

None of the ingredients are listed.

TSCA (Toxic Substances Control Act) :

All ingredients are listed.

CERCLA (Comprehensive Environmental Response, Compensation, and Liability Act):

None of the ingredients are listed.

Proposition 65 (California):

Chemicals known to cause cancer:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients are listed.

Chemicals known to cause developmental toxicity:

None of the ingredients are listed.

Canada

Canadian Domestic Substances List (DSL) :

All ingredients are listed.

SECTION 16: Other information

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.05.2014

Sodium Chloride, ACS Grade

This product has been classified in accordance with hazard criteria of the Controlled Products Regulations and the SDS contains all the information required by the Controlled Products Regulations. Note. The responsibility to provide a safe workplace remains with the user. The user should consider the health hazards and safety information contained herein as a guide and should take those precautions required in an individual operation to instruct employees and develop work practice procedures for a safe work environment. The information contained herein is, to the best of our knowledge and belief, accurate. However, since the conditions of handling and use are beyond our control, we make no guarantee of results, and assume no liability for damages incurred by the use of this material. It is the responsibility of the user to comply with all applicable laws and regulations applicable to this material.

NFPA: 1-0-1

HMIS: 1-0-1

GHS Full Text Phrases: None

Abbreviations and Acronyms:

IMDG International Maritime Code for Dangerous Goods.
PNEC. Predicted No-Effect Concentration (REACH).
CFR Code of Federal Regulations (USA).
SARA Superfund Amendments and Reauthorization Act (USA).
RCRA. Resource Conservation and Recovery Act (USA).
TSCA. Toxic Substances Control Act (USA).
NPRI National Pollutant Release Inventory (Canada).
DOT US Department of Transportation.
IATA International Air Transport Association.
GHS Globally Harmonized System of Classification and Labelling of Chemicals.
ACGIH American Conference of Governmental Industrial Hygienists.
CAS Chemical Abstracts Service (division of the American Chemical Society).
NFPA National Fire Protection Association (USA).
HMIS Hazardous Materials Identification System (USA).
WHMIS Workplace Hazardous Materials Information System (Canada).
DNEL Derived No-Effect Level (REACH).

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

SECTION 1: Identification of the substance/mixture and of the supplier

Product name: Cupric Sulfate, Pentahydrate

Manufacturer/Supplier Trade name:

Manufacturer/Supplier Article number: KEMCU1001-12G

Recommended uses of the product and restrictions on use:

Manufacturer Details:

AquaPhoenix Scientific
9 Barnhart Drive, Hanover, PA 17331
(717) 632-1291

Supplier Details:

AquaPhoenix Scientific Inc.
9 Barnhart Drive, Hanover PA 17331
(717) 632-1291

Emergency telephone number:

Emergency Telephone No.: 800-255-3924

SECTION 2: Hazards identification

Classification of the substance or mixture:



Irritant

Skin irritation, category 2
Eye irritation, category 2A
Acute toxicity (oral), category 4



Environmentally Damaging

Acute hazards to the aquatic environment, category 1
Chronic hazards to the aquatic environment, category 1

Acute Toxicity 4 (oral).
Skin Irritation 2.
Eye Irritation 2.
Aquatic Acute Toxicity 1.
Aquatic Chronic Toxicity 1.

Signal word: Warning

Hazard statements:

Harmful if swallowed.
Causes skin irritation.
Causes serious eye irritation.
Very toxic to aquatic life.
Very toxic to aquatic life with long lasting effects.

Precautionary statements:

If medical advice is needed have product container or label at hand.
Keep out of reach of children.
Read label before use.
Wash thoroughly after handling.
Do not eat, drink or smoke when using this product.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

Wear protective gloves/protective clothing/eye protection/face protection.
Avoid release to the environment.
Do not eat, drink or smoke when using this product.
Rinse mouth.
Specific treatment (see ... on this label).
Take off contaminated clothing and wash before reuse.
Collect spillage.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
IF ON SKIN: Wash with soap and water.
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do - continue rinsing.
If skin irritation or a rash occurs: Get medical advice/attention.
If eye irritation persists get medical advice/attention.

Other Non-GHS Classification: None

SECTION 3: Composition/information on ingredients

Ingredients:

Ingredients:		
CAS 7758-99-8	Copper (II) Sulfate Pentahydrate	100 %
Percentages are by weight		

SECTION 4: First aid measures

Description of first aid measures

After inhalation:

Move exposed individual to fresh air. Loosen clothing as necessary and position individual in a comfortable position. If breathing is difficult, give oxygen. Seek immediate medical attention or advice. Give artificial respiration if necessary.

After skin contact:

Rinse thoroughly. Flush with water for 15 minutes. Seek immediate medical attention or advice.

After eye contact:

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Protect unexposed eye. Rinse/flush exposed eye(s) gently using water for 15-20 minutes. Remove contact lens(es) if able to do so during rinsing. Immediately seek medical assistance or advice.

After swallowing:

Rinse mouth thoroughly. Induce vomiting. Have exposed individual drink sips of water or milk. Seek immediate medical attention or advice.

Most important symptoms and effects, both acute and delayed:

May cause burning of the throat, mouth, esophagus and stomach, nausea, abdominal pain, metallic taste, hemorrhagic gastritis, diarrhea. May cause irritation to the upper respiratory tract, coughing, shortness of breath, sore throat, ulceration and perforation of the respiratory tract, chills and stuffiness of the head. Irritation. Nausea. Headache. Shortness of breath. Copper poisoning may occur if vomiting does not occur immediately. Symptoms may include capillary damage, cold sweat, headache, weak pulse, liver and kidney damage, blood effects, convulsions, jaundice, paralysis, coma, shock or renal failure causing death.

Indication of any immediate medical attention and special treatment needed:

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

Individuals with Wilson's disease are more susceptible to chronic copper poisoning. The use of d-Penicillamine as a chelating agent should be determined by qualified medical personnel. If seeking medical attention, provide SDS document to physician. Note to physician: Treat symptomatically.

SECTION 5: Firefighting measures

Extinguishing media

Suitable extinguishing agents:

If in laboratory setting, follow laboratory fire suppression procedures. Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition.

Unsuitable extinguishing agents: None

Special hazards arising from the substance or mixture:

Combustion products may include metallic oxides or other toxic vapors.

Advice for firefighters:

Protective equipment:

Wear protective equipment. Use respiratory protective device against the effects of fumes, dust, and aerosol. Keep unnecessary personnel away. Ensure adequate ventilation. Keep away from heat and sources of ignition.

Additional information (precautions):

Move product containers away from fire or keep cool with water spray as a protective measure, where feasible.

SECTION 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures:

Wear protective equipment. Avoid dust generation. Use respiratory protective device against the effects of fumes/dust/aerosol. Keep unprotected persons away. Ensure adequate ventilation. Keep away from ignition sources. Protect from heat. Stop the spill, if possible. Contain spilled material by diking or using inert absorbent. Transfer to a disposal or recovery container.

Environmental precautions:

Prevent from reaching drains, sewer or waterway. Collect contaminated soil for characterization per Section 13. Should not be released into environment.

Methods and material for containment and cleaning up:

If in a laboratory setting, follow Chemical Hygiene Plan procedures. Collect liquids using vacuum or by use of absorbents. Place into properly labeled containers for recovery or disposal. If necessary, use trained response staff/contractor.

Reference to other sections: None

SECTION 7: Handling and storage

Precautions for safe handling:

Wash thoroughly after handling. Follow good hygiene procedures when handling chemical materials. Do not eat, drink, smoke, or use personal products when handling chemical substances. If in a laboratory setting, follow Chemical Hygiene Plan. Use only in well ventilated areas. Avoid generation of dust or fine particulate. Wear personal protective equipment. Do not ingest or inhale. Remove contaminated clothing and wash before reuse. Use only in a well-ventilated area. Minimize dust generation and accumulation. Avoid contact with eyes, skin, and clothing. Keep container tightly closed. Avoid ingestion and inhalation.

Conditions for safe storage, including any incompatibilities:

Store away from oxidizing agents. Store in cool, dry conditions in well sealed containers. Store away from incompatible substances. Keep containers tightly closed in a dry, cool and well-ventilated place. Do not expose to air. Store protected from moisture. Store under an inert atmosphere. Store in a cool location. Provide ventilation for containers. Avoid storage near extreme heat, ignition sources or open flame. Store away from foodstuffs.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

SECTION 8: Exposure controls/personal protection



Control parameters:	7758-99-8, Copper sulfate, OSHA PEL TWA: (as Cu) 1 mg/m ³ . 7758-99-8, Copper sulfate, ACGIH TLV: TWA (as Cu) 1 mg/m ³ . 7758-99-8, Copper sulfate, NIOSH REL: TWA (as Cu) 1 mg/m ³ .
Appropriate engineering controls:	Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use/handling. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of vapor or dusts (total/respirable) below the applicable workplace exposure limits (Occupational Exposure Limits-OELs) indicated above.
Respiratory protection:	Not required under normal conditions of use. Use suitable respiratory protective device when high concentrations are present. Use suitable respiratory protective device when aerosol or mist is formed. For spills, respiratory protection may be advisable.
Protection of skin:	The glove material has to be impermeable and resistant to the product/ the substance/ the preparation being used/handled. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.
Eye protection:	Safety glasses with side shields or goggles.
General hygienic measures:	The usual precautionary measures are to be adhered to when handling chemicals. Keep away from food, beverages and feed sources. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Do not inhale gases/fumes/dust/mist/vapor/aerosols. Avoid contact with the eyes and skin.

SECTION 9: Physical and chemical properties

Appearance (physical state, color):	Blue solid	Explosion limit lower:	Not determined
		Explosion limit upper:	Not determined
Odor:	Odorless	Vapor pressure at 20°C:	7.3 mm Hg @ 25 C
Odor threshold:	Not determined	Vapor density:	Not determined
pH-value:	Not determined	Relative density:	Not determined
Melting/Freezing point:	110 C	Solubilities:	Soluble in water.
Boiling point/Boiling range:	150 C	Partition coefficient (n-octanol/water):	Not determined
Flash point (closed cup):	Not applicable	Auto/Self-ignition temperature:	Not determined
Evaporation rate:	Not determined	Decomposition temperature:	Not determined
Flammability (solid, gaseous):	Not determined	Viscosity:	a. Kinematic: Not determined b. Dynamic: Not determined
Density at 20°C:	Not determined		
Specific Gravity:	2.2840 g/cm ³		

SECTION 10: Stability and reactivity

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

Reactivity:

Material is hygroscopic.

Chemical stability:

No decomposition if used and stored according to specifications. Hygroscopic.

Possible hazardous reactions:

None under normal processing.

Conditions to avoid:

High temperatures, dust generation, exposure to moist air or water. Store away from oxidizing agents, strong acids or bases. Incompatible Materials.

Incompatible materials:

Aqueous solution of copper (2+) sulfate is an acid. Strong acids. Strong bases. Incompatible with strong bases, hydroxylamine, magnesium.

Hazardous decomposition products:

Oxides of sulfur. Copper fumes. Carbon oxides (CO, CO₂).

SECTION 11: Toxicological information

Acute Toxicity: No additional information.

Chronic Toxicity: No additional information.

Skin corrosion/irritation:

May cause burning of the throat, mouth, esophagus and stomach.

May cause irritation, redness, pain.

May cause irritation to the upper respiratory tract, ulceration and perforation of the respiratory tract.

Serious eye damage/irritation:

May cause irritation, redness, pain, tearing.

Respiratory or skin sensitization: No additional information.

Carcinogenicity: No additional information.

Germ cell mutagenicity: No additional information.

Reproductive Toxicity: No additional information.

STOT-single and repeated exposure: No additional information.

Additional toxicological information: No additional information.

SECTION 12: Ecological information

Ecotoxicity:

Fish: NOEC (30d) *P. fluviatilis* (7758-98-7), 39 µg/l.

Fish: LC50 (96h) Unspecified Goldfish (7758-98-7), 0.1-2.5 mg/L.

Fish: LC50 (96h) Unspecified Harlequin fish (7758-98-7), 0.1-2.5 mg/L.

Fish: LC50 (96h) Rainbow trout (7758-98-7), 0.1-2.5 mg/L.

Crustacea: EC50 (48h) Unspecified flea *Daphnia* (7758-98-7), 0.24 mg/L.

Fish: LC50 (96 hr) *C. caprio* (7758-98-7 copper sulphate) , 810 µg/l.

Crustacea - LC50 values (48 h) for *D. magna* (7758-98-7), 33.8 µg/L.

Algae: NOEC (72h) *Phaeodactylum tricornutum* (7758-98-7), 5.7 µg/L.

Crustacea: NOEC (27d) *D. polymorpha* (7758-98-7), 21 µg Cu/l.

Persistence and degradability:

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

Readily degradable in the environment.

Bioaccumulative potential:

Copper is accumulated by plants and animals, but it does not appear to biomagnify from plants to animals. This lack of biomagnification appears common with heavy metals. In air, copper aerosols (in general) have a residence time of 2 to 10 days in an unpolluted atmosphere and 0.1 to > 4 days in polluted, urban areas.

Mobility in soil: No additional information.

Other adverse effects: No additional information.

SECTION 13: Disposal considerations

Waste disposal recommendations:

Product/containers must not be disposed together with household garbage. Do not allow product to reach sewage system or open water. It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities (US 40CFR262.11). Consult federal state/ provincial and local regulations regarding the proper disposal of waste material that may incorporate some amount of this product.

SECTION 14: Transport information

US DOT

UN Number:

ADR, ADN, DOT, IMDG, IATA

3288

Limited Quantity Exception:

None

Bulk:

RQ (if applicable): None

Proper shipping Name: Toxic solid, Inorganic, N.O.S., (Cupric Sulfate).

Hazard Class: 6

Packing Group: III.

Marine Pollutant (if applicable): No additional information.

Comments: None

Non Bulk:

RQ (if applicable): None

Proper shipping Name: Toxic solid, Inorganic, N.O.S., (Cupric Sulfate).

Hazard Class: 6

Packing Group: III.

Marine Pollutant (if applicable): No additional information.

Comments: None



SECTION 15: Regulatory information

United States (USA)

SARA Section 311/312 (Specific toxic chemical listings):

Acute

SARA Section 313 (Specific toxic chemical listings):

7758-99-8 Cupric sulfate, Copper compounds.

RCRA (hazardous waste code):

None of the ingredients are listed.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

TSCA (Toxic Substances Control Act) :

All ingredients are listed.

CERCLA (Comprehensive Environmental Response, Compensation, and Liability Act):

7758-99-8 Cupric Sulfate Pentahydrate (Copper compounds) 10 lbs.

Proposition 65 (California):

Chemicals known to cause cancer:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients are listed.

Chemicals known to cause developmental toxicity:

None of the ingredients are listed.

Canada

Canadian Domestic Substances List (DSL) :

All ingredients are listed.

SECTION 16: Other information

For metals listed under CERCLA (antimony, arsenic, beryllium, cadmium, chromium, copper, lead, nickel, selenium, silver, thallium, and zinc), no reporting of releases of the solid form is required if the mean diameter of the pieces of the solid metal released is greater than 100 micrometers (0.004 inches) (Ref: Footnote after Table 302.4 in 40 CFR 302.4). The RQs shown on the consolidated list apply to smaller particles. This product has been classified in accordance with hazard criteria of the Controlled Products Regulations and the SDS contains all the information required by the Controlled Products Regulations. Note. The responsibility to provide a safe workplace remains with the user. The user should consider the health hazards and safety information contained herein as a guide and should take those precautions required in an individual operation to instruct employees and develop work practice procedures for a safe work environment. The information contained herein is, to the best of our knowledge and belief, accurate. However, since the conditions of handling and use are beyond our control, we make no guarantee of results, and assume no liability for damages incurred by the use of this material. It is the responsibility of the user to comply with all applicable laws and regulations applicable to this material.

NFPA: 2-0-1

HMIS: 2-0-1

GHS Full Text Phrases: None

Abbreviations and Acronyms:

IMDG International Maritime Code for Dangerous Goods.

PNEC. Predicted No-Effect Concentration (REACH).

CFR Code of Federal Regulations (USA).

SARA Superfund Amendments and Reauthorization Act (USA).

RCRA. Resource Conservation and Recovery Act (USA).

TSCA. Toxic Substances Control Act (USA).

NPRI National Pollutant Release Inventory (Canada).

DOT US Department of Transportation.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.04.2014

Cupric Sulfate, Pentahydrate

IATA International Air Transport Association.

GHS Globally Harmonized System of Classification and Labelling of Chemicals.

ACGIH American Conference of Governmental Industrial Hygienists.

CAS Chemical Abstracts Service (division of the American Chemical Society).

NFPA National Fire Protection Association (USA).

HMIS Hazardous Materials Identification System (USA).

WHMIS Workplace Hazardous Materials Information System (Canada).

DNEL Derived No-Effect Level (REACH).

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.26.2014

Copper Carbonate

SECTION 1: Identification of the substance/mixture and of the supplier

Product name: Copper Carbonate

Manufacturer/Supplier Trade name:

Manufacturer/Supplier Article number: KEMCC1930-12G

Recommended uses of the product and restrictions on use:

Manufacturer Details:

AquaPhoenix Scientific, Inc.
9 Barnhart Drive
Hanover, PA 17331
1-717-632-1291

Supplier Details:

AquaPhoenix Scientific Inc.
9 Barnhart Drive, Hanover PA 17331
(717) 632-1291

Emergency telephone number:

ChemTel: (24-hour) (US and Canada)
1-(800)-255-3924

SECTION 2: Hazards identification

Classification of the substance or mixture:



Irritant

Acute Tox., oral, 4.
Skin Irrit. 2.
Eye Irrit. 2.
STOT SE 3.
Hazards Not Otherwise Classified - Combustible Dust.

Signal word: Warning

Hazard statements:

Harmful if swallowed.
Causes skin irritation.
Causes serious eye irritation.
May cause respiratory irritation.

Precautionary statements:

If medical advice is needed have product container or label at hand.
Keep out of reach of children.
Read label before use.
Wash skin thoroughly after handling.
Do not eat, drink or smoke when using this product.
Wear protective gloves/protective clothing/eye protection/face protection.
Specific treatment (see supplemental first aid instructions on this label).
Rinse mouth.
Take off contaminated clothing and wash before reuse.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.26.2014

Copper Carbonate

IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
 IF ON SKIN: Wash with soap and water.
 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
 If skin irritation occurs: Get medical advice/attention.
 If eye irritation persists get medical advice/attention.
 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do.
 Continue rinsing.
 Store in a well ventilated place. Keep container tightly closed.
 Dispose of contents and container as instructed in Section 13.

Other Non-GHS Classification: None

SECTION 3: Composition/information on ingredients

Ingredients:

Ingredients:		
CAS 12069-69-1	Copper Carbonate	100 %
Percentages are by weight		

SECTION 4: First aid measures

Description of first aid measures

After inhalation:

Move exposed individual to fresh air. Loosen clothing as necessary and position individual in a comfortable position. Seek medical advice if discomfort or irritation persists. If breathing difficult, give oxygen. If not breathing, give artificial respiration.

After skin contact:

Wash affected area with soap and water. Rinse/flush exposed skin gently using water for 15-20 minutes. Seek medical advice if discomfort or irritation persists. Remove contaminated clothing and shoes. Wash clothing before reuse.

After eye contact:

Protect unexposed eye. Remove contact lens(es) if able to do so during rinsing. Rinse/flush exposed eye(s) gently using water for 15-20 minutes occasionally lifting the upper and lower eyelids. Seek medical assistance.

After swallowing:

Rinse mouth thoroughly. Do not induce vomiting. If victim is conscious and alert, give 2-4 cupfuls of milk or water. Seek medical assistance.

Most important symptoms and effects, both acute and delayed:

Irritation. Nausea. Headache. Shortness of breath. Chronic: Individuals with Wilson's disease are unable to metabolize copper. Thus, copper accumulates in various tissues and may result in liver, kidney, and brain damage.. May cause liver and kidney damage. Aggravated Medical Conditions: Wilson's disease.

Indication of any immediate medical attention and special treatment needed:

If seeking medical attention, provide SDS document to physician. Antidote: The use of d-Penicillamine as a chelating agent should be determined by qualified medical personnel.

SECTION 5: Firefighting measures

Extinguishing media

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.26.2014

Copper Carbonate

Suitable extinguishing agents:

If in laboratory setting, follow laboratory fire suppression procedures. Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition.

Unsuitable extinguishing agents: None

Special hazards arising from the substance or mixture:

Combustion products may include carbon oxides or other toxic vapors. Thermal decomposition can lead to release of irritating gases and vapors. Avoid generating dust; fine dust dispersed in air in sufficient concentrations, and in the presence of an ignition source is a potential dust explosion hazard.

Advice for firefighters:

Protective equipment:

Use NIOSH-approved respiratory protection/breathing apparatus.

Additional information (precautions):

Move product containers away from fire or keep cool with water spray as a protective measure, where feasible. Use spark-proof tools and explosion-proof equipment.

SECTION 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures:

Wear protective equipment. Transfer to a disposal or recovery container. Avoid contact with eyes, skin, and clothing. Use spark-proof tools and explosion-proof equipment. Use respiratory protective device against the effects of fumes/dust/aerosol. Keep unprotected persons away. Ensure adequate ventilation. Keep away from ignition sources. Protect from heat. Stop the spill, if possible. Contain spilled material by diking or using inert absorbent.

Environmental precautions:

Not relevant under normal circumstances considering the small amounts used.

Methods and material for containment and cleaning up:

If in a laboratory setting, follow Chemical Hygiene Plan procedures. Place into properly labeled containers for recovery or disposal. If necessary, use trained response staff/contractor. Dust deposits should not be allowed to accumulate on surfaces, as these may form an explosive mixture if they are released into the atmosphere in sufficient concentration. Avoid dispersal of dust in the air. Collect solids in powder form using vacuum with HEPA filter.

Reference to other sections: None

SECTION 7: Handling and storage

Precautions for safe handling:

Minimize dust generation and accumulation. Wash hands after handling. Avoid dispersal of dust in the air. Routine housekeeping should be instituted to ensure that dusts do not accumulate on surfaces. Dry powders can build static electricity charges when subjected to the friction of transfer and mixing operations. Follow good hygiene procedures when handling chemical materials. Do not eat, drink, smoke, or use personal products when handling chemical substances. Use only in well ventilated areas. Avoid contact with eyes, skin, and clothing.

Conditions for safe storage, including any incompatibilities:

Store in a cool location. Provide ventilation for containers. Avoid storage near extreme heat, ignition sources or open flame. Store away from foodstuffs. Store away from oxidizing agents. Store in cool, dry conditions in well sealed containers. Keep container tightly closed. Store with like hazards.

SECTION 8: Exposure controls/personal protection

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.26.2014

Copper Carbonate



Control parameters:	, , ACGIH TLV TWA (inhalable particles) 10 mg/m ³ . 12069-69-1, Copper Carbonate, OSHA recommended limit 1 mg/m ³ .
Appropriate engineering controls:	Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use/handling. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of vapor or dusts (total/respirable) below the applicable workplace exposure limits (Occupational Exposure Limits-OELs) indicated above. Use under a fume hood. It is recommended that all dust control equipment such as local exhaust ventilation and material transport systems involved in handling of this product contain explosion relief vents or an explosion suppression system or an oxygen deficient environment. Ensure that dust-handling systems (such as exhaust ducts, dust collectors, vessels, and processing equipment) are designed in a manner to prevent the escape of dust into the work area (i.e., there is no leakage from the equipment).
Respiratory protection:	Not required under normal conditions of use. Use suitable respiratory protective device when high concentrations are present. Use suitable respiratory protective device when aerosol or mist is formed. For spills, respiratory protection may be advisable.
Protection of skin:	The glove material has to be impermeable and resistant to the product/ the substance/ the preparation being used/handled. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.
Eye protection:	Safety glasses with side shields or goggles.
General hygienic measures:	The usual precautionary measures are to be adhered to when handling chemicals. Keep away from food, beverages and feed sources. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Do not inhale gases/fumes/dust/mist/vapor/aerosols. Avoid contact with the eyes and skin.

SECTION 9: Physical and chemical properties

Appearance (physical state, color):	Green-blue solid	Explosion limit lower: Explosion limit upper:	Not determined Not determined
Odor:	Odorless	Vapor pressure at 20°C:	Not determined
Odor threshold:	Not determined	Vapor density:	Not determined
pH-value:	Not determined	Relative density:	4.0
Melting/Freezing point:	200 C	Solubilities:	ND = Not Determined.; N/A = Not Applicable.; Insoluble in water.
Boiling point/Boiling range:	24 C	Partition coefficient (n-octanol/water):	Not determined
Flash point (closed cup):	Not determined	Auto/Self-ignition temperature:	Not determined
Evaporation rate:	Not determined	Decomposition temperature:	200 C
Flammability (solid, gaseous):	Not determined	Viscosity:	a. Kinematic: Not determined b. Dynamic: Not determined

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.26.2014

Copper Carbonate

Density at 20°C:	Not determined
-------------------------	----------------

SECTION 10: Stability and reactivity

Reactivity:

Nonreactive under normal conditions.

Chemical stability:

No decomposition if used and stored according to specifications.

Possible hazardous reactions:

None under normal processing. Copper salts + hydrazine react explosively with nitro-methane.

Conditions to avoid:

Store away from oxidizing agents, strong acids or bases. Dust formation, excess heat, incompatible materials.

Incompatible materials:

Strong acids. Strong bases. Strong oxidizers. hydrazine. Copper salts + hydrazine react explosively with nitro-methane.

Hazardous decomposition products:

Copper oxides and carbon oxides. acrid smoke and fumes.

SECTION 11: Toxicological information

Acute Toxicity:

Dermal:

LD50 > 2000 mg/kg bw.

Chronic Toxicity: No additional information.

Skin corrosion/irritation: No additional information.

Serious eye damage/irritation: No additional information.

Respiratory or skin sensitization: No additional information.

Carcinogenicity: No additional information.

Germ cell mutagenicity: No additional information.

Reproductive Toxicity: No additional information.

STOT-single and repeated exposure: No additional information.

Additional toxicological information: No additional information.

SECTION 12: Ecological information

Ecotoxicity: No additional information.

Persistence and degradability:

Readily degradable in the environment.

Bioaccumulative potential: No additional information.

Mobility in soil: No additional information.

Other adverse effects: No additional information.

SECTION 13: Disposal considerations

Waste disposal recommendations:

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.26.2014

Copper Carbonate

Product/containers must not be disposed together with household garbage. Do not allow product to reach sewage system or open water. It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities (US 40CFR262.11). Consult federal state/ provincial and local regulations regarding the proper disposal of waste material that may incorporate some amount of this product.

SECTION 14: Transport information

US DOT

UN Number:

ADR, ADN, DOT, IMDG, IATA

Not Dangerous Goods

Limited Quantity Exception:

None

Bulk:

RQ (if applicable): None

Proper shipping Name: Not Dangerous Goods.

Hazard Class: None

Packing Group: Not Dangerous Goods.

Marine Pollutant (if applicable): No additional information.

Comments: None

Non Bulk:

RQ (if applicable): None

Proper shipping Name: Not Dangerous Goods.

Hazard Class: None

Packing Group: Not Dangerous Goods.

Marine Pollutant (if applicable): No additional information.

Comments: None

SECTION 15: Regulatory information

United States (USA)

SARA Section 311/312 (Specific toxic chemical listings):

Acute

SARA Section 313 (Specific toxic chemical listings):

12069-69-1 Copper Carbonate.

RCRA (hazardous waste code):

None of the ingredients are listed.

TSCA (Toxic Substances Control Act) :

All ingredients are listed.

CERCLA (Comprehensive Environmental Response, Compensation, and Liability Act):

None of the ingredients are listed.

Proposition 65 (California):

Chemicals known to cause cancer:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients are listed.

Safety Data Sheet
according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.26.2014

Copper Carbonate

Chemicals known to cause developmental toxicity:

None of the ingredients are listed.

Canada

Canadian Domestic Substances List (DSL) :

All ingredients are listed.

SECTION 16: Other information

This product has been classified in accordance with hazard criteria of the Controlled Products Regulations and the SDS contains all the information required by the Controlled Products Regulations. Note. The responsibility to provide a safe workplace remains with the user. The user should consider the health hazards and safety information contained herein as a guide and should take those precautions required in an individual operation to instruct employees and develop work practice procedures for a safe work environment. The information contained herein is, to the best of our knowledge and belief, accurate. However, since the conditions of handling and use are beyond our control, we make no guarantee of results, and assume no liability for damages incurred by the use of this material. It is the responsibility of the user to comply with all applicable laws and regulations applicable to this material.

NFPA: 1-0-0

HMIS: 1-0-0

GHS Full Text Phrases: None

Abbreviations and Acronyms:

IMDG International Maritime Code for Dangerous Goods.
PNEC. Predicted No-Effect Concentration (REACH).
CFR Code of Federal Regulations (USA).
SARA Superfund Amendments and Reauthorization Act (USA).
RCRA. Resource Conservation and Recovery Act (USA).
TSCA. Toxic Substances Control Act (USA).
NPRI National Pollutant Release Inventory (Canada).
DOT US Department of Transportation.
IATA International Air Transport Association.
GHS Globally Harmonized System of Classification and Labelling of Chemicals.
ACGIH American Conference of Governmental Industrial Hygienists.
CAS Chemical Abstracts Service (division of the American Chemical Society).
NFPA National Fire Protection Association (USA).
HMIS Hazardous Materials Identification System (USA).
WHMIS Workplace Hazardous Materials Information System (Canada).
DNEL Derived No-Effect Level (REACH).