

# VWR® Swabs for critical environments

01. MULTIPLE HEAD MATERIAL  
OPTIONS FOR A VARIETY OF TASKS

02. NONSTERILE SWABS

03. MADE IN THE USA



**VWR® CRITICAL SWAB® DOUBLE TIPPED COTTON HEAD SWABS**

- Lint-free and residue-free
- Highly absorbent cotton
- Packaged in antistatic leaf pack; 50 leaves of 25 swabs per leaf

An excellent choice when looking for a swab for cleaning sensitive components.



Description	Handle Length, mm (in.)	Handle Width, mm (in.)	Head Length, mm (in.)	Head Width, mm (in.)	Cat. No.	Case of
Cone-Shaped Cotton, Paper Handle	82.0 (3.228)	2.5 (0.100)	11.0 (0.433)	4.7 (0.185)	89031-288	1,250
Slim Cotton, Fine Paper Handle	75.5 (2.972)	1.5 (0.059)	14.0 (0.551)	3.2 (0.126)	89031-292	1,250
Smooth Cotton, Paper Handle	78.0 (3.071)	2.5 (0.100)	15.0 (0.591)	4.4 (0.173)	89031-290	2,500

**VWR® CRITICAL SWAB® COTTON HEAD SWABS**

- Absorbent cotton tip and wood shaft
- Autoclavable
- Disposable

Ideal for general purposes when static or particulate release is less of a concern.



Description	Handle Length, mm (in.)	Handle Width, mm (in.)	Head Length, mm (in.)	Head Width, mm (in.)	Cat. No.	Pack of	Case of
Cotton Head, Wooden Handle	152.0 (6.000)	2.5 (0.100)	15.9 (0.624)	4.8 (0.188)	89031-270	1,000	10,000
Extra Large Cotton Tipped, Wooden Handle	152.0 (6.000)	2.5 (0.100)	17.4 (0.687)	5.9 (0.234)	89031-272	500	5,000

**VWR® CRITICAL SWAB® RETICULATED FOAM SWABS**

- For use in critical environments
- Blue glass-filled polypropylene shaft
- Extra rigid

Thermal edge seals and thermal bonding to the handle ensures that the sponge-like foam cleans up unwanted particles.



Description	Handle Length, mm (in.)	Handle Width, mm (in.)	Head Length, mm (in.)	Head Width, mm (in.)	Cat. No.	Pack of	Case of
Elongated Foam Tipped, Polypropylene Handle	82.8 (3.260)	2.5 (0.100)	25.4 (1.000)	3.6 (0.140)	89031-284	500	1,000
Cone-Shaped Reticulated Foam Head, Polypropylene Handle	83.0 (3.270)	2.5 (0.100)	17.5 (0.688)	3.9 (0.154)	89031-282	500	1,000
Fine-Point Reticulated Foam Head, Polypropylene Handle	161.0 (6.324)	2.4 (0.094)	14.6 (0.574)	1.8 (0.070)	89031-294	500	1,000

**VWR® CRITICAL SWAB® LARGE RECTANGULAR FOAM HEAD SWAB**

- 100% polyurethane foam head with 100ppi
- Thermally bonded foam head reduces contamination
- Sturdy polypropylene handle



Description	Handle Length, mm (in.)	Handle Width, mm (in.)	Head Length, mm (in.)	Head Width, mm (in.)	Cat. No.	Pack of	Case of
Rectangular Foam Tip, Polypropylene Handle	130.0 (5.100)	5 (0.195)	25.4 (1.000)	13.5 (0.531)	89031-286	500	1,000

**VWR® CRITICAL SWAB® MEDIUM FOAM HEAD SWABS**

- 100% polyurethane foam head with 100ppi
- Wood shaft
- Disposable



89031-274

Description	Handle Length, mm (in.)	Handle Width, mm (in.)	Head Length, mm (in.)	Head Width, mm (in.)	Cat. No.	Pack of	Case of
Medium Foam Head, Wooden Handle	152.0 (5.984)	2.5 (0.100)	17.4 (0.685)	4.8 (0.188)	89031-274	500	1,000

**VWR® CRITICAL SWAB® PADDLE-SHAPED FOAM HEAD SWAB**

- 100% polyurethane foam head with 100ppi
- Blue polypropylene shaft
- Flexible and disposable



89031-276

Description	Handle Length, mm (in.)	Handle Width, mm (in.)	Head Length, mm (in.)	Head Width, mm (in.)	Cat. No.	Pack of	Case of
Paddle-Shapped Foam Head, Polypropylene Handle	152.0 (5.984)	3.1 (0.120)	15.1 (0.594)	4.8 (0.188)	89031-276	500	1,000

**VWR® CRITICAL SWAB® FOAM-OVER-COTTON HEAD SWAB**

- Foam-over-cotton tips made from shred-resistant polyurethane foam over USP pharma-grade cotton
- Great for cleaning with solvents
- Wood shaft
- Disposable



89031-280

Description	Handle Length, mm (in.)	Handle Width, mm (in.)	Head Length, mm (in.)	Head Width, mm (in.)	Cat. No.	Pack of	Case of
Foam-Over-Cotton, Wooden Handle	152.0 (6.000)	2.5 (0.100)	23.8 (0.937)	9.5 (0.375)	89031-280	500	1,000

**VWR® CRITICAL SWAB® MINI FOAM HEAD SWAB**

- Thermal edge seals and thermal bonding to the handle
- Sponge-like foam cleans up unwanted particles without releasing any into your critical environment
- Black polypropylene handle



89176-450

Description	Handle Length, mm (in.)	Overall Length, mm (in.)	Handle Diameter, mm (in.)	Tip Diameter, mm (in.)	Tip Size, mm (in.)	Cat. No.	Pack of	Case of
Foam Tip, Polypropylene Handle	66.6 (2.625)	68.58 (2.700)	2.28 (0.09)	3.96 (0.156)	15.23 (0.6)	89176-450	500	2,500

# Setting science in motion to create a better world

At Avantor®, a new era has begun. From breakthrough discovery to agile delivery of advanced products and services across the globe, Avantor manufactures and distributes high-quality products, services and solutions to professionals in the life sciences and advanced technologies industries.



VWR, part of Avantor provides an integrated, seamless purchasing experience optimized for your success.

Learn how the new Avantor is moving forward at [www.settingscienceinmotion.com](http://www.settingscienceinmotion.com)

©2019 Avantor, Inc. All rights reserved.



800 932 5000 | [VWR.COM](http://VWR.COM)

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.