


















































































































































Every glove is not the same. Below, we've recommended gloves that match your needs

- Alternative Energy
- Biotechnology
- Biopharmaceutical
- Contract Research/Contract Manufacturing
- Clinical/Healthcare
- Educational and Medical Research
- Environmental Testing
- Food and Beverage
- Government
- Histology
- Med Device
- Microelectronics
- Nanotechnology
- Petrochemical
- Pharmaceutical
- Water/Wastewater
- University

	Microflex® Product	Microflex® Item #	VWR Item #	Size	Gloves Per Box	Description	Alternative Energy	Biotechnology	Biopharmaceutical	Contract Research/Contract Manufacturing	Clinical/Healthcare	Educational and Medical Research	Environmental Testing	Food and Beverage	Government	Histology	Med Device	Microelectronics	Nanotechnology	Petrochemical	Pharmaceutical	Water/Wastewater	University
NATURAL RUBBER LATEX	 ¹	CTP-233	89107-530	XS-L	100	High quality polymer coated latex, fully textured, comfortable latex formula, \$1 per case donated to BCRF																	
	 ¹	EV-2050	32916-534	XS-XL	100	Premium latex glove, maximum comfort, dons easily in wet and dry environments, gold standard for Microflex latex products																	
	 ¹	MF-300	32916-500	XS-XL	100	Latex glove, thicker at the fingertips, benchmark in strength and durability																	
	 ¹	DGP-350	32916-636	XS-XL	100	All purpose latex glove, comfortable fit, good tactile sensitivity, fully textured for a secure grip																	
	 ¹	SG-375	32916-602	S-XL	50	Extra thick extended cuff latex gloves are the best choice in a latex for heavy-duty/high-risk applications																	
	 ¹	CFG-900	14900-536	XS-XL	100	Best overall value in a quality latex glove, fully textured for a secure grip																	
POWDER-FREE	 ²	DFK-608	89094-570	S-2XL	50	Nitrile glove, ultra-thin moisture absorbing flock-lining, formulated to reduce hand fatigue, fully-textured for a firm grip																	
	 ²	TQ-601	94024-060	XS-XL	100	Tranquility™ gloves are soft, thin, and comfortable for everyday use. Tranquility™ gloves are latex-free and contain no natural rubber proteins. Compact box size saves on storage requirements.																	
	 ²	SU-690	32916-662	XS-XL	100	Nitrile glove, best combination of strength and durability in an all purpose synthetic.																	
	 ²	SEC-375	32916-672	S-3XL	50	Nitrile glove, best combination of strength and durability in an all purpose synthetic. Tested for use with chemotherapy drugs.																	
	 ²	FFS-700	32916-684	XS-XL	100	Nitrile gloves are the best overall value in a comfortable mid-range synthetic glove																	
	 ²	US-220	32919-004	XS-XL	100	Nitrile glove with an extended cuff, patented formulation provides comfort and protection																	
	 ²	MK-296	94001-368	XS-3XL	100	Black nitrile glove, fully textured, distinct color masks stains																	
	 ²	NPG-888	32916-792	XS-XL	100	Chloroprene gloves are micro-textured for working in wet environments and are tested for use with chemotherapy drugs																	
VINYL *	 ²	DF-850	32916-804	XS-XL	100	Powder-free vinyl glove, latex-free formulation with a smooth finish.																	
































* Microflex® Nitrile/Synthetic gloves are not made from natural rubber proteins.

¹ Caution (LATEX): This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



² Caution (SYNTHETIC): Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

1 U.S. Patent No. 6,451,893, U.S. Patent No. 7,176,260

Synthetic and Clean Room

Synthetic and Clean Room						Alternative Energy	Biotechnology	Biopharmaceutical	Contract Research/ Contract Manuf.	Clinical/Healthcare	Educational and Medical Research	Environmental Testing	Food and Beverage	Government	Histology	Med Device	Microelectronics	Nanotechnology	Petrochemical	Pharmaceutical	Water/Wastewater	University	
Microflex® Product	Microflex® Item #	VWR Item #	Size	Gloves Per Box	Description																		
NATURAL RUBBER LATEX		DC-400	32916-618	XS-XL	1000	Lightly powdered vinyl glove, latex-free formulation with a smooth finish																	
		CE5-429	32916-704	S-XL	1000	9" powder-free industrial latex cleanroom glove, class 100																	
		CE5-512	32916-712	XS-XL	1000	12" powder-free industrial latex cleanroom glove, class 100. Low surface particulates, ionic extractables and non-volatile residues. Micro-textured fingertip surface enhances wet or dry grip. Silicone free, provides superb dexterity, tactile sensitivity, fit, feel and comfort.																	
		CE4-200	32916-720	XS-XL	1000	12" powder-free industrial latex cleanroom glove, class 10. Washed in microprocessor controlled, two-in-one, washer/dryer laundering system to provide very low surface particulates, ionic extractables and non-volatile residues. Micro textured surface enhances wet or dry grip. Silicone free, provides superb dexterity, tactile sensitivity, fit, feel and comfort. Double vacuum sealed packaging.																	
NITRILE * 1		CE5-637	32916-756	XS-XL	1000	9" powder-free industrial nitrile cleanroom glove, class 100. Low surface particulates, ionic extractables and non-volatile residues. Silicone free, provides superb dexterity, tactile sensitivity, fit, feel and comfort.																	
		CE5-755	32916-764	XS-XL	1000	9" powder-free industrial nitrile cleanroom glove, class 100. Low surface particulates, ionic extractables and non-volatile residues. Silicone free, provides superb dexterity, tactile sensitivity, fit, feel and comfort.																	
		CE4-313	32916-772	XS-XL	1000	12" powder-free industrial nitrile cleanroom glove, class 10. Washed in microprocessor controlled, two-in-one, washer/dryer laundering system to provide very low surface particulates, ionic extractables and non-volatile residues. Does not contain plasticizers. Silicone free, provides superb dexterity, tactile sensitivity, fit, feel and comfort. Excellent static decay and surface resistivity properties. Double vacuum sealed packaging.																	
HAND SPECIFIC LATEX		HSCE4-879	32916-746	6.0-9.0	1000	12" powder-free industrial latex cleanroom glove, class 10. Washed in microprocessor controlled, two-in-one, washer/dryer laundering system to provide very low surface particulates, ionic extractables and non-volatile residues. Micro textured surface enhances wet or dry grip. Silicone free, provides superb dexterity, tactile sensitivity, fit, feel and comfort. Hand specific design minimizes unnecessary stress and strain that can cause musculoskeletal injury over time. Irradiation indicators are placed on all master poly-bags and cartons - your guarantee of sterility assurance																	

* Microflex® Nitrile/Synthetic gloves are not made from natural rubber proteins.
 1 U.S. Patent No. 6,451,893,
 U.S. Patent No. 7,176,260

 **Caution (LATEX):** This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.
 **Caution (SYNTHETIC):** Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



1.800.932.5000 vwr.com

Prices, product appearance and specifications are current at the time of printing, subject to change without notice. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in U.S. dollars unless otherwise noted. Offers valid in USA and Canada, void where prohibited by law or company policy, while supplies last. Visit vwr.com to view our privacy policy and additional disclaimers.
 VWR, forms of VWR, and the VWR logo and/or design are either registered trademarks ®, trademarks ™, or service marks ™ of VWR International, LLC in the United States and/or other countries. All other marks referenced herein are registered trademarks, trademarks or service marks of their respective owner(s). For a complete list of trademark owners please visit vwr.com.
 ©2010 VWR International, LLC. All rights reserved.