

Think. Inside. The. Box.

PURELAB® flex 2*



DI Pack

LC208

Lifetime expectancy: 12 months
What is purified from the water?
Ions

Impact of not changing:
Weakly bound ions are released back into the water reducing resistivity



UV light (185nm)

LC210

Lifetime expectancy: 12-18 months
What is purified from the water?
Purifies TOC, inactivates DNA and RNA
Impact of not changing:

Uncontrolled bacterial contamination and TOC levels are not reduced



Point of Use Filter**

LC145

Lifetime expectancy: 3 months
What is purified from the water?
Particulates and bacteria
Impact of not changing:

Potential recontamination of purified water and reduced flow rates

*A combination of treatments and technologies produce the final water specification

** For DNase, RNase & endotoxin specification at point of use, please use a LC197 Biofilter