



Item Number C1939
Item Calcium Gluconate, Anhydrous, Powder, USP
CAS Number 299-28-5
Molecular Formula $C_{12}H_{22}CaO_{14}$
Molecular Weight 430.37

Assay.....98.0-102.0%

MAXIMUM LIMITS

Loss on Drying..... 3.0%
Chloride..... 0.07%
Sulfate.....0.05%
Arsenic.....3 ppm
Heavy Metals.....0.002%
Reducing Substances..... 1.0%
Residual Solvents.....To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: 14422 S. San Pedro St.
Gardena, CA 90248
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843

East Coast Facility: 755 Jersey Ave.
New Brunswick, NJ 08901