

Axygen® Axypet® Pro Pipettors

Color-coded for your comfort



Axygen Axypet Pro Pipettors have been engineered to provide the highest levels of comfort, accuracy, and precision. The lightweight construction, contoured handle, and 4-digit counter were designed to ensure comfortable pipetting. All pipettors feature smooth plunger movement and extremely low pipetting forces to reduce wrist strain and fatigue (RSI). Colored pushbuttons enable easy size identification and reduce the risk of cross-contamination.

- ▶ Ergonomic handle design is a perfect fit for both the right and left hand
- ▶ Convenient one-hand volume setting with auto-lock system prevents accidental volume change during work
- ▶ 4-digit counter provides greater precision and is perfectly positioned for visibility at all times
- ▶ Single-channel pipettors are compatible with very narrow tubes
- ▶ Fully autoclavable and UV-resistant for sample protection against contaminations
- ▶ The multi-channel shafts retract individually for perfect tip loading and easier tip ejection
- ▶ Unique ejector lever system enables effortless tip ejection
- ▶ Easy in-lab calibration helps you quickly adjust your pipettor to various liquid types
- ▶ Individually tested and supplied with a certificate of quality

Axygen Axypet Pro Starter Kit

Axygen Axypet Pro Starter Kit contains four variable volume pipettors in the most popular and universal volumes (0.5-10 μ L, 2-20 μ L, 20-200 μ L, and 100-1000 μ L), pipettor stand for four pipettors, and 3 racks with tips.



Ordering Information

Axygen® Axypet® Pro Pipettors

Single-channel

VWR Cat. No.	Corning Cat. No.	Volume Range (µL)	Accuracy (%)	Precision (%)	Non-filtered Tip (µL)	Qty/Cs
76320-080	AP-2-P	0.1 - 2	±40.0 to ±1.5	≤12.0 to ≤0.7	10	1
76320-082	AP-10-P	0.5 - 10	±4.0 to ±1.0	≤2.8 to ≤0.4	10	1
76320-084	AP-20-P	2 - 20	±3.5 to ±0.8	≤1.5 to ≤0.3	200	1
76320-086	AP-50-P	5 - 50	±3.5 to ±0.8	≤1.3 to ≤0.3	200	1
76320-088	AP-100-P	10 - 100	±3.0 to ±0.8	≤1.0 to ≤0.2	200	1
76320-090	AP-200-P	20 - 200	±2.0 to ±0.6	≤0.7 to ≤0.2	200	1
76320-092	AP-1000-P	100 - 1000	±2.5 to ±0.6	≤0.6 to ≤0.2	1000	1

8-channel

76320-094	AP-8-10-P	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
76320-096	AP-8-50-P	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
76320-098	AP-8-200-P	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
76320-100	AP-8-300-P	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1

12-channel

76320-102	AP-12-10-P	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
76320-104	AP-12-50-P	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
76320-106	AP-12-200-P	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
76320-108	AP-12-300-P	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1

Axygen Axypet Pro Starter Kit

VWR Cat. No.	Corning Cat. No.	Description	Qty/Cs
76320-110	AP-STR-KIT-P	Axygen Axypet Pro Starter Kit Includes: 4 single-channel pipettors: 0.5 to 10 µL (Cat. No. 76320-082), 2 to 20 µL (Cat. No. 76320-084), 20 to 200 µL (Cat. No. 76320-090), 100 to 1000 µL (Cat. No. 76320-092); universal linear stand for four single-channel pipettors (Cat. No. 76324-848); Axxygen MultiRack pipet tips: 10 µL (Cat. No. 76320-676), 200 µL (Cat. No. 76320-678), 1000 µL (Cat. No. 76320-684)	1

Accessories

VWR Cat. No.	Corning Cat. No.	Description	Qty/Cs
76324-848	AP-USD4-P	Universal linear stand for four pipettors, red (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
76324-852	AP-USD6	Universal linear stand for six pipettors, transparent (six single-channel, maximum of three 8-channel or 12-channel pipettors)	1

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.


800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.