

TECHNIGLOVE

INTERNATIONAL



ORDERING INFORMATION

TechnipakTM 100-count plastic boxes
CLASS 100 COMPATIBLE
 9-1/2" powder-free blue nitrile gloves

		Cat. No.
TNI00PFB	X-SMALL	89415-066
TNI01PFB	SMALL	89415-068
TNI02PFB	MEDIUM	89415-070
TNI03PFB	LARGE	89415-072
TNI04PFB	X-LARGE	89415-074
TNI05PFB	XX-LARGE	89415-076

100 Gloves/TechnipakTM, 10 TechnipaksTM/Case = 1,000 Gloves/Case, Case Size: 11"x13"x11"

TechniToteTM 250-count plastic boxes
CLASS 100 COMPATIBLE
 12" blue nitrile gloves

		Cat. No.
TNTI200PFB	X-SMALL	89415-054
TNTI201PFB	SMALL	89415-056
TNTI202PFB	MEDIUM	89415-058
TNTI203PFB	LARGE	89415-060
TNTI204PFB	X-LARGE	89415-062
TNTI205PFB	XX-LARGE	89415-064

250 Gloves/TechniToteTM, 4 TechniTotesTM/Case = 1,000 Gloves/Case, Case Size: 14"x16"x13"

PHYSICAL PROPERTIES

Style: Non-Sterile	Ambidextrous
Lengths: Technipak TM	9-1/2" (240mm)
TechniTote TM	12" (300mm)
Thickness:	5 mil (.005")
Grip Surface:	Microtextured Fingertips
Cuff:	Beaded
Color:	Blue
Tensile Strength:	>12.5 MPa
Elongation at break: (%)	>500

TECHNICAL PROPERTIES

Particle Levels: Tote/Pak	<2500 total particles/cm ² >0.5µm IEST - RP-CC005.2 Method
ESD Properties:	At 50-60% Humidity
Surface Resistivity:	<10-10 ohm/square per ESD-S11:11
Static Decay:	<5 seconds per RETS-5-003

Extractables

TNVR:	<5.0 µg/cm ²	Phosphate:	<0.05 µg/cm ²
Chloride:	<1.5 µg/cm ²	Fluoride:	<0.05 µg/cm ²
Sulfate:	<1.0 µg/cm ²	Nitrate:	<0.05 µg/cm ²
Nitrate:	<1.0 µg/cm ²	Bromide:	<0.05 µg/cm ²