

**FEATURES**

The molded seat provides comfort as well as durability and easy cleaning on this **standard height chair**. Fully adjustable flip-up padded arms provide comfort and safety.



**Blood Draw Chair**

Item# 12CPS

**ACCESSORIES**




Drawer/Tray Combo



Flip-up Tray




Clear Arm Wrap



Optional Straight Arms

**SPECIFICATIONS:**

- Contoured Seat
- 6" Adjustable Arm Height
- Extra-wide Padded "L" Arms



Seat Height	Patient Width	Weight Capacity	Overall Dimensions
18"	20"	300 lbs	34"W x 24"D

**READY TO SHIP PRODUCT: 12CPS**

This chair is currently not available as a Ready To Ship item.

**CERTIFICATION(S)**

As an option we offer CAL-133 protection, a flammability test procedure for seating furniture for use in Public Occupancies that has been adopted by many states. The test results show that the product passed both Criteria Q and Criteria B for the California TB-133 Burn Test.

