



Product Specification

Catalog Number: BDH9298

Catalog Description: BDH SOD PHOSPHAT MONO ACS

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics

Specifications

Calcium	<= 0.005 %
Chloride	<= 0.0005 %
Heavy Metals	<= 0.001 %
Insolubles	<= 0.01 %
Iron	<= 0.001 %
pH (5%, Water) @25C	4.1 - 4.5
Potassium	<= 0.01 %
Purity	98.0 - 102.0 %
Sulfate	<= 0.003 %

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager