

Certificate of Quality



Part Number	MFLX06435-24
Lot Number	29206023
Part Description	TUBING MFLEX PHARMAPUR NO. 24 25FT
Expiration Date	01FEB2025
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: PharmaPure
- Meets USP Class VI requirements.
- Registered with a Drug Master File with the FDA.
- These products have been manufactured in an ISO-9001 registered facility.
- A lot control system is used to assure traceability to all raw materials.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.