

Safety Data Sheet

According to Hazardous Products Regulation (SOR/2015-17)

Revision date: 20.03.2024

Version: 1.0

Print date: 20.03.2024

SECTION 1: Identification

Product identifier

Trade name/designation:	IC Instrument Ion Anion Check Standard 2 in H2O
Product No.:	BDH0833-125ML
Synonyms:	none
CAS No.:	Not applicable
Other means of identification:	

Relevant identified uses of the substance or mixture and uses advised against

Recommended use:	For Further Manufacturing Use Only
Uses advised against:	Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Supplier

VWR International

Street	2360 Argentia Road
Postal code/City	Mississauga, Ontario Canada L5N 5Z7
Telephone	+1-800-932-5000 toll-free within US/Canada
Telefax	+1-610-728-2103

Emergency phone number

Telephone	+1-613-996-6666 (Canutec, 24 hrs/day, 7 days/week, Canada)
-----------	------------------------------------------------------------

Preparation Information

VWR International - Product Information Compliance

E-mail SDS@avantorsciences.com

SECTION 2: Hazard identification

2.1 Classification of the substance or mixture

Classification according to Hazardous Products Regulation (SOR/2015-17)

This mixture is classified as not hazardous according to Hazardous Products Regulation (SOR/2015-17)

2.2 Label elements

Labelling in accordance with (SOR/2015-17)

According to Hazardous Products Regulation (SOR/2015-17) the product does not have to be labelled.

Hazard(s) not otherwise classified (HNOC)
none

SECTION 3: Composition/information on ingredients

3.1 Substances

not applicable

3.2 Mixtures

Hazardous ingredients GHS Classification in accordance with (SOR/2015-17)

Substance name	Concentration	Identifier	Hazard classes and hazard categories
Sodium nitrate	< 0.1%	CAS No.: 7631-99-4	Ox. Sol. 3 - H272 Eye Irrit. 2 - H319
Sodium fluoride	< 0.1%	CAS No.: 7681-49-4	Acute Tox. 3 - H301 Skin Irrit. 2 - H315 Eye Irrit. 2 - H319

SECTION 4: First aid measures

4.1 Description of first aid measures

General information

When in doubt or if symptoms are observed, get medical advice. Change contaminated, saturated clothing. Wash contaminated clothing before reuse. Do not leave affected person unattended.

In case of inhalation

Remove casualty to fresh air and keep warm and at rest. Obtain medical attention if symptoms appear.

In case of skin contact

Gently wash with plenty of soap and water. In case of skin reactions, consult a physician.

After eye contact:

Rinse immediately carefully and thoroughly with eye-bath or water. Obtain medical attention if symptoms appear.

In case of ingestion

Rinse mouth thoroughly with water. Call a doctor if you feel unwell.

Self-protection of the first aider

First aider: Pay attention to self-protection!

4.2 Most important symptoms/effects, acute and delayed

No known symptoms to date.

4.3 Indication of any immediate medical attention and special treatment needed

No special information on medical attention and special treatment available.

SECTION 5: Fire fighting measures

5.1 Extinguishing media

Suitable extinguishing media

The product itself does not burn.

Co-ordinate fire-fighting measures to the fire surroundings.

Extinguishing media which must not be used for safety reasons

Full water jet.

5.2 Specific hazards arising from the chemical

In case of fire may be liberated:

Pyrolysis products, toxic

5.3 Advice for firefighters

The product itself does not burn.

Do not breathe smoke.

Protective equipment and precautions for firefighters:

Wear a self-contained breathing apparatus and chemical protective clothing.

Do not allow run-off from fire-fighting to enter drains or water courses.

Do not inhale explosion and combustion gases.

Use water spray jet to protect personnel and to cool endangered containers.

In case of fire: Evacuate area.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

For non-emergency personnel: Remove victim out of the danger area. First Aid, decontamination, treatment of symptoms.

6.2 Environmental precautions

No special environmental measures are necessary.

6.3 Methods and material for containment and cleaning up

Take up mechanically, placing in appropriate containers for disposal. Dispose according to legislation.

6.4 Reference to other sections

Personal protection equipment (PPE): see section 8 Disposal information: see section 13

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advices on safe handling

No special measures are necessary.

Measures to prevent fire, aerosol and dust generation

No special measures are necessary.

Measures required to protect the environment

No special measures are necessary.

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

7.2 Conditions for safe storage, including any incompatibilities

Recommended storage temperature: no data available

Store in a well-ventilated place. Keep container tightly closed. Packaging materials: High density polyethylene (HDPE) Unsuitable container/equipment material: Metal container

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Does not contain substances above concentration limits fixing an occupational exposure limit.

8.2 Engineering controls

Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment. If handled uncovered, arrangements with local exhaust ventilation have to be used.

Personal protection equipment (PPE)

Wear suitable protective clothing. When handling with chemical substances, protective clothing must be worn.

Eye/face protection

Eye glasses with side protection

Skin protection

Wear suitable gloves. When handling with chemical substances, protective gloves must be worn. In the case of wanting to use the gloves again, clean them before taking off and air them well. Check leak tightness/impermeability prior to use.

By short-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,12 mm
Breakthrough time	> 480 min

By long-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,12 mm
Breakthrough time	> 480 min

Respiratory protection

Respiratory protection necessary at: aerosol or mist formation If exposure limits are exceeded or irritation is experienced, NIOSH approved respiratory protection should be worn.

Additional information

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

Environmental exposure controls

no data available

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance	
Physical state:	liquid
Color:	colorless
Odor:	characteristic

Safety relevant basic data

pH:	no data available
Melting point/freezing point:	no data available
Initial boiling point and boiling range:	100 °C (1013 hPa)
Flash point:	no data available
Flammability:	Not applicable
Lower and upper explosion limit	
Lower explosion limit:	no data available
Upper explosion limit:	no data available
Vapor pressure:	23 hPa (20 °C)
Relative vapour density:	no data available
Density and/or relative density	
Density:	1.00118 g/cm ³ (20 °C)
Solubility(ies)	
Water solubility:	soluble (20°C)
Partition coefficient: n-octanol/water:	no data available
Auto-ignition temperature:	no data available
Decomposition temperature:	Not applicable
Viscosity	
Kinematic viscosity:	no data available
Dynamic viscosity:	no data available
Particle characteristics:	does not apply to liquids

9.2 Other information

Evaporation rate:	no data available
Explosive properties:	no data available
Oxidising properties:	Not applicable
Bulk density:	no data available
Refraction index:	no data available
Dissociation constant:	no data available
Surface tension:	no data available
Henry's Law Constant:	no data available

SECTION 10: Stability and reactivity

10.1 Reactivity

This material is non-reactive under normal conditions.

10.2 Chemical stability

The product is chemically stable under standard ambient conditions (room temperature).

10.3 Possibility of hazardous reactions

No information available.

10.4 Conditions to avoid

Keep away from sources of heat (e.g. hot surfaces), sparks and open flames.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous decomposition products

No known hazardous decomposition products.

Decomposition products in case of fire: see section 5.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute effects

Acute oral toxicity:

Sodium nitrate - LD50: > 1267 mg/kg - Rat - (RTECS)

Sodium fluoride - LD50: > 31 mg/kg - Rat - (RTECS)

Sodium fluoride - LDLo: > 71 mg/kg - Human - (RTECS)

Acute dermal toxicity:

Sodium fluoride - LD50: 175 mg/kg - Rat - (National Library of Medicine ChemID Plus (NLM CIP))

Acute inhalation toxicity:

no data available

Irritant and corrosive effects:

Primary irritation to the skin:

Not applicable

Irritation to eyes:

Not applicable

Irritation to respiratory tract:

Not applicable

Respiratory or skin sensitization

In case of skin contact: not sensitizing

In case of inhalation: not sensitizing

STOT-single exposure

Not applicable

STOT-repeated exposure

Not applicable

CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**Carcinogenicity**

No indication of human carcinogenicity.

Germ cell mutagenicity

No indications of human germ cell mutagenicity exist.

Reproductive toxicity

No indications of human reproductive toxicity exist.

Aspiration hazard

Not applicable

Other adverse effects

no data available

Additional information

no data available

SECTION 12: Ecological information

12.1 Toxicity

Fish toxicity:

Sodium nitrate - LC50: 573 - 10000 mg/l (96 h) - McGurk, M.D., F. Landry, A. Tang, and C.C. Hanks 2006. Acute and Chronic Toxicity of Nitrate to Early Life Stages of Lake Trout (*Salvelinus namaycush*) and Lake Whitefish (*Coregonus clupeaformis*). *Environ.Toxicol.Chem.* 25(8):2187-2196

Sodium fluoride - LC50: 51 - 830 mg/l (96 h) - Saxena, R., R. Gupta, M. Tripathi, and K. Gopal 2001. Fluoride Induced Haematological Alterations in the Freshwater Fish *Channa punctatus*. *J.Ecophysiol.Occup.Health* 1(1/2):139-146

Daphnia toxicity:

Sodium nitrate - LC50: 3400 - 5600 mg/l (48 h) - Dowden, B.F., and H.J. Bennett 1965. Toxicity of Selected Chemicals to Certain Animals. *J.Water Pollut.Control Fed.* 37(9):1308-1316

Sodium fluoride - EC50: 98 mg/l (48 h) - Dave, G. 1984. Effects of Fluoride on Growth, Reproduction and Survival in *Daphnia magna*. *Comp.Biochem.Physiol.C* 78(2):425-431

Sodium fluoride - LC50: 340 mg/l (48 h) - LeBlanc, G.A. 1980. Acute Toxicity of Priority Pollutants to Water Flea (*Daphnia magna*). *Bull.Environ.Contam.Toxicol.* 24(5):684-691 (OECDG Data File)

Algae toxicity:

Sodium fluoride - EC50: 900 mg/l (96 h) - Kuehn, R., and M. Pattard 1990. Results of the Harmful Effects of Water Pollutants to Green Algae (*Scenedesmus subspicatus*) in the Cell Multiplication Inhibition Test. Water Res. 24(1):31-38 (OECDG Data File)

Sodium fluoride - EC50: 850 mg/l (72 h) - Kuehn, R., and M. Pattard 1990. Results of the Harmful Effects of Water Pollutants to Green Algae (*Scenedesmus subspicatus*) in the Cell Multiplication Inhibition Test. Water Res. 24(1):31-38 (OECDG Data File)

Bacteria toxicity:

no data available

12.2 Persistence and degradability

no data available

12.3 Bioaccumulative potential

Partition coefficient: n-octanol/water: no data available

12.4 Mobility in soil:

no data available

12.5 Results of PBT/vPvB assessment

Not applicable

12.6 Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to the environment.

12.7 Other adverse effects

no data available

SECTION 13: Disposal considerations**13.1 Waste treatment methods****Appropriate disposal / Product**

Dispose according to legislation. Consult the appropriate local waste disposal expert about waste disposal.

Waste code product: no data available

Appropriate disposal / Package

Dispose according to legislation. Handle contaminated packages in the same way as the substance itself.

Additional information

none

No further relevant information available.

SECTION 14: Transport information

Land transport (TDG)

14.1	UN number or ID number:	No dangerous good in sense of this transport regulation.
14.2	UN proper shipping name:	not assigned
14.3	Transport hazard class(es):	none
14.4	Packing group:	not assigned
14.5	Environmental hazards:	none
14.6	Special precautions for user:	none

Sea transport (IMDG)

14.1	UN number or ID number:	No dangerous good in sense of this transport regulation.
14.2	UN proper shipping name:	not assigned
14.3	Transport hazard class(es):	none
14.4	Packing group:	not assigned
14.5	Environmental hazards:	none
14.6	Special precautions for user:	none
14.7	Maritime transport in bulk according to IMO instruments	not relevant

Air transport (ICAO-TI / IATA-DGR)

14.1	UN number or ID number:	No dangerous good in sense of this transport regulation.
14.2	UN proper shipping name:	not assigned
14.3	Transport hazard class(es):	none
14.4	Packing group:	not assigned
14.5	Special precautions for user:	none

SECTION 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

Domestic Substance List:

SECTION 16: Other information

Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists
DOT - Department of Transportation
IARC - International Agency for Research on Cancer
IATA-DGR - International Air Transport Association-Dangerous Goods Regulations
ICAO-TI - International Civil Aviation Organization-Technical Instructions
IMDG - International Maritime Code for Dangerous Goods
LTV - Long Term Value
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety & Health Administration
PBT - Persistent, Bioaccumulative and Toxic
PEL - Permissible Exposure Limit
STV - Short Term Value
SVHC - Substances of Very High Concern
TDG - Transport of Dangerous Goods
TLV - Threshold Limit Value
vPvB - very Persistent, very Bioaccumulative

Key literature references and sources for data

This Safety Data Sheet has been prepared based on information available for public as TOXNET information, European Chemicals Agency (ECHA) substance dossier, papers from international cancer research institutes (IARC Monographs), U.S. National Toxicology Program data, U.S. Agency for Toxic Substances and Disease Control (ATSDR), PubChem websites and SDS from our raw material manufacturers.

Revision date	Version	Print date
20.03.2024	1.0	20.03.2024

Additional information

Indication of changes none/none

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guidance. The information in this document is based on the present state knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and his Affiliates shall not be held liable for any damage resulting from handling.