



WHEATON®

CRIMP/SERUM VIALS



DWK Life Sciences is the manufacturer of high-quality medical grade glass vials formed from the highest quality pharmaceutical-grade tubing, conforming to USP <660> / EP3.2.1 / JP <7.01> standards for injectable packaging and the safe storage and delivery of drug content. Vials are designed for the short- and long-term storage of parenteral, lyophilized, reagent, and vaccine content.

Product Highlights

- | | |
|---|--|
| <ul style="list-style-type: none">• Clear and amber, Type I, tubular, low-expansion, borosilicate glass conforms to ASTM and USP requirements | <ul style="list-style-type: none">• Global ISO 15378 certified manufacturing |
| <ul style="list-style-type: none">• ISO 8362-1, GPI 2710, and Custom design | <ul style="list-style-type: none">• Certified pre-sterilized Ready-to-Use configurations |
| <ul style="list-style-type: none">• European, American, and no blowback* | <ul style="list-style-type: none">• CCI validated stoppers, seals, and snap caps |
| <ul style="list-style-type: none">• Small minimum order requirements | <ul style="list-style-type: none">• Autoclavable |

*No blowback is only available in a custom made vial

DWK Life Sciences SALES PARTNER

GPI 2710 EP US Blowback Serum Vial Configurations

Catalog Pack
Shrink Wrapped Tray

Carton Bulk Pack
Shrink-wrapped Trays in Case Packs

Unitized Pallet
Shrink-wrapped Trays on Unitized Pallets

VWR Standard Pack Cat. No.	Size (mL)	Brimful Capacity (mL)	Catalog Pack Qty	VWR Bulk Pack Cat. No.	Carton	Case Qty	VWR Pallet Cat. No.	Vials/ Unitized Pallet
Clear								
76716-344	2	3.2	225	77392-006	12	2,700	77392-020	68,400
76716-346	3	4.8	169	77392-008	10	1,690	77392-022	51,376
76716-348	5	10	100	77392-012	10	1,000	77392-026	30,400
76716-350	10	15.2	150	77392-016	14	2,100	77392-030	39,600
Amber								
76716-336	2	3.2	225	77391-990	12	2,700	77391-998	68,400
76716-338	3	4.8	169	77391-992	10	1,690	77392-000	51,376
76716-340	5	10	100	77391-994	10	1,000	77392-002	30,400
76716-342	10	15.2	150	77391-996	14	2,100	77392-004	39,600

GPI 2710 EP US Blowback Serum Vial Dimensions

VWR Cat. No.	OD (mm)	Height (mm)	Wall (mm)	Finish (mm)	Tray Size (inches)	Carton	Carton Size (inches)	Pallet
2 mL Clear								
77392-020	14.75	32	1.00	13	9.125 x 7.75 x 1.25	—	—	Plastic
77392-006	14.75	32	1.00	13	9.125 x 7.75 x 1.25	12	16.125 x 9.875 x 8.5625	—
76716-344	14.75	32	1.00	13	9.125 x 7.75 x 1.25	1	9.875 x 8.125 x 2.125	—
2 mL Amber								
77391-998	14.75	32	1.00	13	9.125 x 7.75 x 1.25	—	—	Plastic
77391-990	14.75	32	1.00	13	9.125 x 7.75 x 1.25	12	16.125 x 9.875 x 8.5625	—
76716-336	14.75	32	1.00	13	9.125 x 7.75 x 1.25	1	9.875 x 8.125 x 2.125	—
3 mL Clear								
77392-022	16.75	37.7	1.10	13	9.125 x 7.75 x 1.25	—	—	Plastic
77392-008	16.75	37.7	1.10	13	9.125 x 7.75 x 1.25	10	16.125 x 9.875 x 8.5625	—
76716-346	16.75	37.7	1.10	13	9.125 x 7.75 x 1.25	1	9.875 x 8.125 x 2.125	—
3 mL Amber								
77392-000	16.75	37.7	1.10	13	9.125 x 7.75 x 1.25	—	—	Plastic
77391-992	16.75	37.7	1.10	13	9.125 x 7.75 x 1.25	10	16.125 x 9.875 x 8.5625	—
76716-338	16.75	37.7	1.10	13	9.125 x 7.75 x 1.25	1	9.875 x 8.125 x 2.125	—
5 mL Clear								
77392-026	21.75	40	1.20	20	9.125 x 7.75 x 1.25	—	—	Plastic
77392-012	21.75	40	1.20	20	9.125 x 7.75 x 1.25	10	16.125 x 9.875 x 8.5625	—
76716-348	21.75	40	1.20	20	9.125 x 7.75 x 1.25	1	9.875 x 8.125 x 2.125	—
5 mL Amber								
77392-002	21.75	40	1.20	20	9.125 x 7.75 x 1.25	—	—	Plastic
77391-994	21.75	40	1.20	20	9.125 x 7.75 x 1.25	10	16.125 x 9.875 x 8.5625	—
76716-340	21.75	40	1.20	20	9.125 x 7.75 x 1.25	1	9.875 x 8.125 x 2.125	—
10 mL Clear								
77392-030	23.75	50	1.30	20	12.5 x 10.5 x 1.75	—	—	Plastic
77392-016	23.75	50	1.30	20	12.5 x 10.5 x 1.75	14	21.625 x 13.25 x 15.0625	—
76716-350	23.75	50	1.30	20	12.5 x 10.5 x 1.75	1	13.25 x 10.875 x 3.3125	—
10 mL Amber								
77392-004	23.75	50	1.30	20	12.5 x 10.5 x 1.75	—	—	Plastic
77391-996	23.75	50	1.30	20	12.5 x 10.5 x 1.75	14	21.625 x 13.25 x 15.0625	—
76716-342	23.75	50	1.30	20	12.5 x 10.5 x 1.75	1	13.25 x 10.875 x 3.3125	—

ISO 8362-1 EP European Blowback Serum Vial

Catalog Pack
Shrink Wrapped Tray

Carton Bulk Pack
Shrink-wrapped Trays in Case Packs


Unitized Pallet
Shrink-wrapped Trays on Unitized Pallets

VWR Standard Pack Cat. No.	Size (mL)	Brimful Capacity (mL)	Catalog Pack Qty	VWR Bulk Pack Cat. No.	Carton	Case Qty	VWR Pallet Cat. No.	Vials/Unitized Pallet
Clear								
77392-188	2	4	196	77392-010	12	2,352	77392-024	59,584
77392-190	6	10	100	77392-014	10	1,000	77392-028	30,400
77392-192	10	13.5	150	77392-018	16	2,400	77392-032	43,200

ISO 8362-1 EP European Blowback Serum Vial

VWR Cat. No.	OD (mm)	Height (mm)	Wall (mm)	Finish (mm)	Tray Size (inches)	Carton	Carton Size (inches)	Pallet
2R Clear Vial								
77392-024	16	35	1.00	13	9.125 x 7.75 x 1.25	—	—	Plastic
77392-010	16	35	1.00	13	9.125 x 7.75 x 1.25	12	16.125 x 9.875 x 8.5625	—
77392-188	16	35	1.00	13	9.125 x 7.75 x 1.25	1	—	—
6R Clear Vial								
77392-028	22	40	1.00	20	9.125 x 7.75 x 1.25	—	—	Plastic
77392-014	22	40	1.00	20	9.125 x 7.75 x 1.25	10	16.125 x 9.875 x 8.5625	—
77392-190	22	40	1.00	20	9.125 x 7.75 x 1.25	1	—	—
10R Clear Vial								
77392-032	24	45	1.00	20	12.5 x 10.5 x 1.75	—	—	Plastic
77392-018	24	45	1.00	20	12.5 x 10.5 x 1.75	16	21.625 x 13.25 x 15.0625	—
77392-192	24	45	1.00	20	12.5 x 10.5 x 1.75	1	—	—





Excellence in your hands

038 / WTD 11-2023

© 2024 DWK Life Sciences • The trademarks used are owned by DWK Life Sciences GmbH, DWK Life Sciences LLC or DWK Life Sciences Ltd.

DWK Life Sciences SALES PARTNER

 **avantor**[™]
delivered by **VWR**[™]

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2024 Avantor, Inc. All rights reserved.

0324 Lit. No. 040090W