



Material	37964-690
Material description	Ether anhydrous
Grade	ACS
CAS Number	60-29-7
Molecular formula	(CH <sub>3</sub> CH <sub>2</sub> ) <sub>2</sub> O
Molecular mass	74.12
Storage	Store between 15° C and 25° C.

## Characteristics

## Specifications

Assay	Min. 99.0 %
Ethanol	Passes test
Stabilized	Passes test
Color (Pt-Co)	Max. 10 APHA
Titration acid	Max. 0.0002 meq/g
Carbonyl (as HCHO)	Max. 0.001 %
Peroxide (as H <sub>2</sub> O <sub>2</sub> )	Max. 1 ppm
Residue after evaporation	Max. 0.001 %
Water	Max. 0.03 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.