



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1127-4LG, BDH1127-4LP, BDH1127-19L, BDH1127-204L, CABDH1127-4LG, CABDH1127-4LP, CABDH1127-19L, CABDH1127-204L,

Catalog Description: 99% n-Heptane

Grade:

Date Current Specification Set: April 19 2006

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

| Characteristics | Specifications |
|-----------------|---|
| Appearance | Clear |
| Water (K.F.): | Less than 0.010% |
| Residue: | Less than 0.001% (0.0010 g/146 mL) |
| GC Purity: | Greater than 99.9% n-heptane (n-C7) and C7 isomers |
| n-Heptane: | Greater than 99.0% |
| Color (APHA): | Less than 10 |

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor