

Product Name:	VWR-HC-FRP-17 76579-150
Description:	<p>VWR Catalog Number: 76579-150  VWR Standard Series Natural Refrigerant Flammable Storage Refrigerator  17 Cu. Ft. capacity  One solid swing door, right hinged, not reversible  Microprocessor Temperature Controller  Adjustable operating temperature range: 1°-10°C  Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms  Sample simulated digital temperature display; .1° Resolution  Adjustable temperature display in Celsius or Fahrenheit  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  5 adjustable shelves  Keyed door lock  Magnetic door gasket for positive seal  Heavy duty door handle  Leveling legs  One probe access port (3/8")  White textured cabinet and door  High density urethane foam cabinet and door insulation  Natural, Hydrocarbon Refrigerant (R600a)  Two year parts and labor warranty  Utilizes ATEX and non-sparking components for safety  Exterior dimensions: 30 3/4" W x 28 5/8"D x 67 3/4"H  Cycle defrost  115V, 60 Hz, 3 Amps, 1/4 HP  ETL listed  Units meet NFPA guidelines 45 and 70  Shipping weight: 240 lbs.</p>