





REVISION NO	REVISION DATE	DESCRIPTION	REVISED BY	APPROVED BY

ROUND WIDE MOUTH BOTTLE SPECIFICATION

VWR-US CAT NO	VWR EUROPE CAT NO	MATERIAL	COLOR	CAPACITY	APROX BRIM CAPACITY (ml)	CAP SIZE	HEIGHT WITH OUT CAP (mm)	BODY DIA (mm)
414004-186	215-1889	HDPE	TRANSLUCENT	4L	4.5L	100-415	298	163
414004-187	215-1898	PP	TRANSLUCENT	2L	2.28L	100-415	240	119
414004-188	215-1899	PP	TRANSLUCENT	4L	4.5L	100-415	298	163
414004-205	215-1918	HDPE	TRANSLUCENT	2L	2.28L	100-415	240	119

NOTE :
1) PP MATERIAL BOTTLES ARE AUTOCLABLE.
2) 4L BOTTLE HAS INTEGRATED HANDLE DESIGN.



PRODUCT IMAGE FOR 2L CAPACITY BOTTLES



PRODUCT IMAGE FOR 4L CAPACITY BOTTLES

TITLE : WIDE MOUTH ROUND BOTTLE

DRAWN BY	
DRAWN DATE	2/2/2009
DRAWING NO	

Material : AS PER TABLE

Rev : 0

SCALE : NOT TO SCALE

