

# Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Cuff Tape
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2021-10-20

VWR Catalogue Number	Dimensions	Sterility
89495-444	2" x 36 yds	Non - Sterile
89495-446	2" x 36 yds	Sterile

## Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Product Specifications

<b>Material</b>	White perforated polyethylene tape.
<b>Sterilization</b>	Products labeled as Sterile is gamma irradiated and dose released upon EN ANSI/AAMI/ISO 11137 recommended practices in effect at the time of validation. Product labeled as sterile meets a minimum and maximum specified dose range of 40 kGy - 60 kGy respectively.
<b>ATP Assay</b>	Not applicable.
<b>DNase &amp; RNase Free</b>	Not applicable.
<b>BSE/TSE</b>	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

<b>Cytotoxicity</b>	Not applicable.
<b>Non-Pyrogenic Statement</b>	Not applicable.
<b>Latex Statement</b>	Not applicable.
<b>BPA Statement</b>	Not applicable.
<b>DEHP Statement</b>	Not applicable.
<b>RoHS</b>	Not applicable.
<b>REACH Statement</b>	Not applicable.
<b>Storage Conditions</b>	Avoid exposure to direct sunlight, UV light, extreme temperature, and humidity fluctuations.
<b>Shelf Life</b>	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.