


AC-1262xx	CUSTOMER DRAWING	Page 1 of 1
Revision 02	Cryo.s™ 2ml, round bottom, starfoot base	 greiner bio-one
	Item-No. 126 2xx	

NOT TO SCALE



Item No.	Description	A	B	C	D	E	F
126 261	Cryo.s™ 2ml	Ø 9.8	Ø 12.4	Ø 13.1	40	47.3	14.5
126 263	Cryo.s™ 2ml	Ø 9.8	Ø 12.4	Ø 13.1	40	47.3	14.5
126 277	Cryo.s™ 2ml	Ø 9.8	Ø 12.4	Ø 13.1	40	47.3	14.5
126 278	Cryo.s™ 2ml	Ø 9.8	Ø 12.4	Ø 13.1	40	47.3	14.5
126 279	Cryo.s™ 2ml	Ø 9.8	Ø 12.4	Ø 13.1	40	47.3	14.5
126 280	Cryo.s™ 2ml	Ø 9.8	Ø 12.4	Ø 13.1	40	47.3	14.5

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date November 20, 2008	Date November 20, 2008	Date November 21, 2008	
Date December 22, 2004	Name F. Falter	Name F. Sonntag	Name F. Sonntag	