

Quality Certificate

VWR® HDPE MULTIPURPOSE CONTAINERS

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	89009-662
Manufacture Lot Number:	43071 20200831
Description:	VWR Container 4OZ W/LID CS250
Date of Manufacture:	August 31, 2020
Country of Origin:	China

Quality System Compliance

Products are manufactured under the ISO 9001:2015 & ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing

Representative production samples are collected and inspected in accordance with current applicable standard protocols.

Dimensional Analysis: ASQ Z1.4-2003, General Insp. Level- II, AQL- 0.65. **Result:** Within tolerance

Product Specifications

This product is Class I category device as defined by the FDA in 21CFR Parts 862-892.

Material:

Container: High Density Polyethylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed

Lid: Blend of High Density Polyethylene and Polyolefin Elastomer, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed

Expiration date: Expiration date (Shelf life) is based on product and packaging shelf life provided the product wrapper is not damaged.

BSE/TSE Statement: All materials used and contact materials are of fully synthetic origin.

BPA Statement: N/A

DEHP Statement: N/A

Non-Pyrogenic: N/A

DNase & RNase Free: N/A

Cytotoxicity: N/A

Latex free Statement: N/A.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 14, 2020