



## VWR® Sterile Polycarbonate Single-Use Flask Assemblies

- Sterile flasks come individually bagged and process ready
- Multiple flask sizes allow for process specific flexibility
- When used with sterile VWR tube sets, custom possibilities are limitless

### Sterile Polycarbonate flasks with weldable tubing provide the maximum amount of flexibility and reliability

- Assembled in ISO Class 7 cleanroom
- Comprehensive size selection
- Tubing is completely weldable/sealable
- USP Class VI materials
- All product contact materials are animal-free or processed in accordance to EMEA/410/01
- Able to withstand temperatures between -20°C and 65°C without losing integrity
- Assemblies individually doubled bagged with product label located on inner bag
- Gamma irradiated to Sterility Assurance Level (SAL) 10<sup>-6</sup>

### Assemblies have been designed with a broad spectrum of applications in mind

- Life Sciences & General Lab Use
- Environmental Assays
- Production & Cleanroom Applications

## VWR® Sterile Polycarbonate Single-Use Flask Assemblies

| Capacity, mL (oz.) | Cap Size | Tubing Material | Tube Size (ID x OD; Length)                                                     | End Connection                                                                | Pack of | Case of | Cat. No.  |
|--------------------|----------|-----------------|---------------------------------------------------------------------------------|-------------------------------------------------------------------------------|---------|---------|-----------|
| 125 (4)            | 38-430   | C-Flex          | Port 1: 1/8" x 1/4"; 24"<br>Port 2: 1/8" x 1/4"; 3"                             | Port 1: PP Tube Plug<br>Port 2: 25mm PES Vent Filter                          | 6 Each  | 4 Packs | 10830-700 |
| 250 (8)            | 38-430   | C-Flex          | Port 1: 1/8" x 1/4"; 24"<br>Port 2: 1/8" x 1/4"; 3"                             | Port 1: PP Tube Plug<br>Port 2: 25mm PES Vent Filter                          | 6 Each  | 4 Packs | 10830-702 |
| 500 (16)           | 38-430   | C-Flex          | Port 1: 1/8" x 1/4"; 24"<br>Port 2: 1/8" x 1/4"; 3"                             | Port 1: PP Tube Plug<br>Port 2: 25mm PES Vent Filter                          | 6 Each  | 4 Packs | 10830-278 |
| 1000 (32)          | 53B      | C-Flex          | Port 1: 3/8" x 1/2"; 12"<br>Port 2: 1/4" x 3/8"; 3"<br>Port 3: 3/8" x 1/2"; 12" | Port 1: PP Tube Plug<br>Port 2: 50mm PTFE Vent Filter<br>Port 3: PP Tube Plug | 3 Each  | 4 Packs | 10830-280 |
| 2000 (64)          | 53B      | C-Flex          | Port 1: 3/8" x 1/2"; 12"<br>Port 2: 1/4" x 3/8"; 3"<br>Port 3: 3/8" x 1/2"; 12" | Port 1: PP Tube Plug<br>Port 2: 50mm PTFE Vent Filter<br>Port 3: PP Tube Plug | 2 Each  | 6 Packs | 10830-290 |

Consult VWR for custom stopper and bottle options; contact your VWR Sales Representative or call 1.800.932.5000.