

# VWR<sup>®</sup> Safety Cabinets for Flammables and Corrosives

- 01. FULLY WELDED, 18-GAUGE DOUBLE-WALL CONSTRUCTION FOR FIRE PROTECTION
- 02. UPPER AND LOWER DIAMETRICALLY OPPOSED SIDE VENTS
- 03. ALL MODELS FIRE TESTED AND APPROVED BY FM GLOBAL



VWR® Safety Cabinets are designed to protect both laboratory personnel and the cabinet contents. Fully welded, 18-gauge steel ensures fire protection and strength with reduced gaps. Equipped with upper and lower diametrically opposed side vents, a ground connector on lower right side of the cabinet, double-wall construction with 3.8cm (1½") insulating air space, and an environmentally friendly, electrostatically applied powder coat finish.

All door styles meet OSHA 1910.106, NFPA 30 and NFPA 1.

Self-close door styles also meet IFC requirements and feature a fusible link that melts at 74°C (165°F) for automatic closure in the event of a fire.



### For Flammables

| D x W x H, cm (in.)                 | Door Type  | No. of Shelves | Cat. No.  |
|-------------------------------------|------------|----------------|-----------|
| <b>30 Gallon Storage</b>            |            |                |           |
| 45.7 x 109.2 x 111.8 (18 x 43 x 44) | Self Close | 1              | 89522-616 |
| 45.7 x 109.2 x 111.8 (18 x 43 x 44) | Bi-Fold    | 1              | 89522-618 |
| 45.7 x 109.2 x 111.8 (18 x 43 x 44) | Manual     | 1              | 89522-620 |
| <b>45 Gallon Storage</b>            |            |                |           |
| 45.7 x 109.2 x 165 (18 x 43 x 65)   | Self Close | 2              | 89522-622 |
| 45.7 x 109.2 x 165 (18 x 43 x 65)   | Bi-Fold    | 2              | 89522-624 |
| 45.7 x 109.2 x 165 (18 x 43 x 65)   | Manual     | 2              | 89522-626 |
| <b>60 Gallon Storage</b>            |            |                |           |
| 86.3 x 86.3 x 165 (34 x 34 x 65)    | Self Close | 2              | 89522-628 |
| 86.3 x 86.3 x 165 (34 x 34 x 65)    | Bi-Fold    | 2              | 89522-630 |
| 86.3 x 86.3 x 165 (34 x 34 x 65)    | Manual     | 2              | 89522-632 |
| <b>90 Gallon Storage</b>            |            |                |           |
| 86.3 x 109.2 x 165 (34 x 43 x 65)   | Self Close | 2              | 89522-634 |
| 86.3 x 109.2 x 165 (34 x 43 x 65)   | Bi-Fold    | 2              | 89522-636 |
| 86.3 x 109.2 x 165 (34 x 43 x 65)   | Manual     | 2              | 89522-638 |

### For Corrosives

| D x W x H, cm (in.)                 | Door Type  | No. of Shelves | Cat. No.  |
|-------------------------------------|------------|----------------|-----------|
| <b>30 Gallon Storage</b>            |            |                |           |
| 45.7 x 109.2 x 111.8 (18 x 43 x 44) | Self Close | 1              | 89522-640 |
| 45.7 x 109.2 x 111.8 (18 x 43 x 44) | Manual     | 1              | 89522-642 |
| <b>45 Gallon Storage</b>            |            |                |           |
| 45.7 x 109.2 x 165 (18 x 43 x 65)   | Self Close | 2              | 89522-644 |
| 45.7 x 109.2 x 165 (18 x 43 x 65)   | Manual     | 2              | 89522-646 |
| <b>60 Gallon Storage</b>            |            |                |           |
| 86.3 x 86.3 x 165 (34 x 34 x 65)    | Self Close | 2              | 89522-648 |
| 86.3 x 86.3 x 165 (34 x 34 x 65)    | Manual     | 2              | 89522-650 |

### Accessories

| Description                                  | Cat. No.  |
|--|-----------|
| Spillslope Shelf for 30 and 45 Gal. Cabinets | 89522-652 |
| Spillslope Shelf for 60 Gal. Cabinet         | 89522-654 |
| Spillslope Shelf for 90 Gal. Cabinet         | 89522-656 |

