


# PurFlock® Ultra Flocked Swab

Standard tip, polystyrene handle with breakpoint at 80mm



## Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US 8420385

Patents pending: Canada 2,744,678, Europe 10251389.2

## Dimensions

**6.04" Overall Length** (153.416 mm)

### Tip

■ **0.2" Dia.** (5.08 mm)

■ **0.69" Length** (17.526 mm)

### Handle

■ **0.1" Dia.** (2.54 mm) tapered to 0.066" (1.676 mm)

■ **6" Length** (152.654 mm)

## Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of sterility and/or conformance available with each shipment upon request

## Qualities

■ Medical grade quality

■ Highly absorbent tip

■ Soft smooth tip

## Applications

■ For collection of bacteria or micro-organisms

■ Buccal cell collection

■ Specimen collection

■ Forensic evidence collection

■ Rapid diagnostic testing

■ DNA testing

Ref #	Packaging
<b>3306-U</b> <i>Non-sterile</i>	10/100 (10 bags per case, 100 per bag)
<b>25-3306-U</b> <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 swab per package)



**PurFlock®  
Ultra**

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

**Tel** 800-321-2313 (US and Canada Only)  
207-876-3311

**Fax** 800-323-4153  
207-876-3130

**Web** [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

**Email** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
**[puritanmedproducts.com](http://puritanmedproducts.com)**