

Product Name:	HCUCBI0420ADLH
Description:	<p> VWR Catalog No. 76210-114 VWR Undercounter Freezer (Built-In)-LEFT HINGED 4.2 Cu. Ft. Capacity One solid swing door Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: -15°C – -25°C Two fixed shelves Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 23 ¾ "W x 24 ½" D x 31 15/16 "H Manual defrost 115V, 60 Hz, 1.5 Amps, 1/8 HP ENERGY STAR CERTIFIED ADA compliant UL/C-UL listed Shipping weight: 133 lbs. </p>