

Product Name:	HCUCBI0404SS
Description:	<p>VWR Catalog No. 10819-886</p> <p>VWR Undercounter Refrigerator (Built-In)</p> <p>4.5 Cu. Ft. Capacity</p> <p>One solid stainless steel swing door, right hinged* Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Audible and visual high and low temperature alarms</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Remote alarm contacts</p> <p>Adjustable operating temperature range: 1°C – 10°C</p> <p>Three shelves (two adjustable/one fixed)with guard rail on back to prevent product from pushing against back wall</p> <p>Back wall evaporator cover</p> <p>Forced draft circulation</p> <p>Recessed door handle</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>Stainless steel exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty</p> <p>Exterior dimensions: 23 ¾ "W x 24 ½" D x 33 3/8 "H</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 1.74 Amps, 1/8 HP</p> <p>ENERGY STAR CERTIFIED</p> <p>UL/C-UL listed</p> <p>Shipping weight: 135 lbs.</p> <p>*This product is available left hinged, product number HCUCBI0404SSLH</p>