

VORTEX-GENIE PULSE Pulsing Vortex Mixer



Vortex-Genie Pulse Accessories

More Than Just Vortexing!

- True pulsing action provides more collisions and random motion.
- For applications requiring more aggressive action than standard vortexing.
- Variable "RUN/STOP" time provides quick pulsing or interval mixing.
- Fully programmable "ON" and "REST" function reduces heat buildup within the vessel, while keeping materials in suspension. Ideal for heat sensitive samples.
- Use with all Vortex-Genie accessories.
- Ideal for applications such as re-suspensions of lyophilized materials, creating solutions, enzymatic reactions and much more.

Specifications

Supplied	3-inch Platform, Pop-off Cup
Speed	500 - 3000 RPM (2850 for 50Hz models)
Timer - Hands Free	1-99 minutes or continuous
Pulse - ON/OFF	0.1 seconds - 99 minutes
Weight	4 kgs (8.8 lbs)
Dimensions (D x W x H)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)



Cat. No.	Description	Cat. No.	Description
89202-298	Vortex-Genie Pulse	58815-180	Rubber Cover only for 3-inch Platform
58815-230	Model H301 Multiple Sample Starter Set	89202-322	TurboMix™ Attachment
58815-212	6-inch Platform	89202-324	Universal Microtube Holder - (Requires 89202-322)
58815-206	Microtube Foam Insert (2) - (Requires 58815-212)	13502-388	Large Ampule/Tube Attachment
58815-185	Microplate Foam Insert (2) - (Requires 58815-212)	13502-390	Small Ampule/Tube Attachment
89415-150	Deepwell Microplate Foam Insert (2) - (Requires 58815-212)	97060-174	Horizontal 50ml Tube Holder (6 tubes)
58815-224	Model H302 Large Sample Set	97060-176	Horizontal 15ml Tube Holder (12 tubes)
58815-222	Recessed Platform without Bands	97060-178	Horizontal Microtube Holder (24 tubes)
58815-188	9-16mm Tube Foam Inserts (2) - (Requires 58815-222)	89401-416	Vertical High Speed 50ml Tube Holder (3 tubes)
58815-186	14-29mm Tube Foam Inserts (2) - (Requires 58815-222)	89401-418	Vertical High Speed 15ml Tube Holder (6 tubes)
58815-184	29-37mm Tube Foam Inserts (2) - (Requires 58815-222)	89401-420	Vertical High Speed Microtube Holder (9 tubes)
13502-386	Blank Foam Inserts (2) - (Requires 58815-222)	97060-180	Vertical 50ml Tube Holder (6 tubes)
58815-220	Assorted Foam Inserts (4) - (Requires 58815-222)	97060-182	Vertical 15ml Tube Holder (12 tubes)
58815-218	Elastic Bands (2) - (Requires 58815-222)	97060-184	Vertical Microtube Holder (24 tubes)
58815-210	Pop-off Cup	97060-186	Vertical (6) 50ml, (6) 15ml, (12) Microtube Holder
14221-972	3-inch Platform with Rubber Cover	89202-338	Bottom Closure with Screws Kit