

VWR[®] SHAKING WATER BATHS

Stability $\pm 0.1^{\circ}\text{C}$

Robust and reliable

User friendly with easy
to operate controls

Amb. $+5^{\circ}\text{C}$ to 99°C
operation

User settable sample
protection

Tray and lid
included as
standard

3-year
warranty





Specifications

		 12L	 18L
Tank capacity	L/gal.	12 / 3.2	18 / 4.8
Shaking speed	rpm	20 - 200	
Linear shaking	stroke length / mm	20	
User calibration		single point	
Over temperature protection		user-settable	
Overall consumption	120V	0.8kW	1.05kW
Ext. dimensions	h x d x w / cm	27.5 x 38 x 36	27.5 x 56.5 x 33.5
	h x d x w / in	10.8" x 14.9" x 14.1"	10.8" x 22.2" x 13.1"
Shaking tray area	cm	24 x 23.5	42 x 23.5
	in	9.4" x 9.3"	16.5" x 9.3"
Lid included		●	
Catalog number		10128-126	10128-128

Accessories (full range available on website)

Test tube tray

Compatible with test tube racks. Choice of 8 racks to accommodate different tube diameters and microtubes.



Description	Catalog no
Test tube tray 12L bath (holds up to 3 test tube racks)	10144-848
Test tube tray 18L bath (holds up to 5 test tube racks)	10144-850

Rack type	Tube size Ø	Capacity	Catalog no
Microtube rack	0.5 ml	119	10144-836
Microtube rack	1.5 ml	48	10144-838
Test tube rack	10 mm	48	10144-828
Test tube rack	13 mm	44	10144-830
Test tube rack	16 mm	24	10144-852
Test tube rack	19 mm	21	10144-854
Test tube rack	25 mm	12	10144-856
Test tube rack	30 mm	10	10144-858

Replacement universal tray with springs

Designed for ultimate versatility as to the types and sizes of vessel which can be accommodated. Adjustable spring configuration for maximum flask capacity.



Description	Catalog no
Universal tray for 12L bath	10144-844
Universal tray for 18L bath	10144-846

Stainless steel base tray

Allows bath to be used as an unstirred bath



Description	Catalog no
Base tray for 12L bath	10144-840
Base tray for 18L bath	10144-842

Note: Baths include non-drip polycarbonate lid and universal tray as standard.



1.800.932.5000 | vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.