


PDS No. 1272xx	PRODUCT DATA SHEET					Page 1 of 1
Revision 05	Cryo.s™, 4 ml, Round Bottom, Starfoot Base				CE	 greiner bio-one
	Item-No. 1272xx					
Valid for Item-No.:	127261 (sterile)	127263 (sterile)	127277 (sterile)	127278 (sterile)	127279 (sterile)	
	127280 (sterile)					

1.	Description / Specification	
1.1	Description	Cryo.s™ with different coloured screw caps, external thread, round bottom, starfoot base and additional cap inserts for labelling 127261: without writing area 127263, -277, -278, -279, -280: with writing area and graduation Additional inserts, Cat.-No. 304 171 (50 pieces per bag)
1.2	Dimensions	13,1 x 83,3 mm (with screw cap)
1.3	Volume	Working volume: up to 4,0 ml
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Screw cap: HDPE (High Density Polyethylene), free of heavy metal Small cap insert: PS (Polystyrene), free of heavy metal
1.5	Colour	Tube: natural, writing area and graduation: white Screw cap: 127261, -263: natural 127277: green 127278: yellow 127279: blue 127280: red Small cap insert: white
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	-196°C to +100°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	5800 x g: swinging-bucket rotor 34000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: <a href="http://www.gbo.com/bioscience">www.gbo.com/bioscience</a> →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	5 years after month of production
2.7	Other Information	Recommended for storage in the freezer or above liquid nitrogen

3.	Packaging	
3.1	Pieces / Bag	50
3.2	Pieces / Box	300
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Additional in each box: Certificate of Quality Safety precaution 2 bags with 50 cap inserts per bag (100 pieces) 1 felt pen

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 16 May 2013	Date 17 May 2013	Date 17 May 2013	
Date 11.12.2009	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name A. Schulz	