

Consistent Quality



VWR® Dual Click Specimen Containers

VWR's Dual Click line of specimen containers are manufactured from virgin medical grade propylene. Containers feature two audible and tactile "clicks" to let the user know the lid is properly secured. Dual Click closures ensure safe transportation of specimen. Once the specimen is collected, a simple turn of the closure will give the user assurance that the seal is fully engaged with a "click" sound. The second "click" ensures 95 kPa compliance for pneumatic tube transportation. Manufactured and sterilized (sterile containers) in the USA.

Volume	Width	Cap Color	Cat. No.		Case of
			Sterile	Non-Sterile	
60 mL	48 mm	Natural	10803-950	10803-998	500
60 mL	48 mm	White	10803-952	10804-000	500
60 mL	48 mm	Green	10803-954	10804-002	500
60 mL	48 mm	Blue	10803-956	10804-004	500
60 mL	48 mm	Yellow	10803-958	10804-006	500
60 mL	48 mm	Orange	10803-960	10804-008	500
60 mL	53 mm	White	10836-000	10836-002	560
90 mL	48 mm	Natural	10803-962	10804-010	400
90 mL	48 mm	White	10803-964	10804-012	400
90 mL	48 mm	Green	10803-966	10804-014	400
90 mL	48 mm	Blue	10803-968	10804-016	400
90 mL	48 mm	Yellow	10803-970	10804-018	400
90 mL	48 mm	Orange	10803-972	10804-020	400
90 mL	53 mm	Natural	10803-974	10804-022	400
90 mL	53 mm	White	10803-976	10804-024	400
90 mL	53 mm	Green	10803-978	10804-026	400
90 mL	53 mm	Blue	10803-980	10804-028	400
90 mL	53 mm	Yellow	10803-982	10804-030	400
90 mL	53 mm	Orange	10803-984	10804-032	400
120 mL	53 mm	Natural	10803-986	10804-034	300
120 mL	53 mm	White	10803-988	10804-036	300
120 mL	53 mm	Green	10803-990	10804-038	300
120 mL	53 mm	Blue	10803-992	10804-040	300
120 mL	53 mm	Yellow	10803-994	10804-042	300
120 mL	53 mm	Orange	10803-996	10804-044	300
Individually Wrapped, Sterile					
90 mL	48 mm	Natural	10804-052	-	100
90 mL	53 mm	Natural	10804-050	-	100
120 mL	53 mm	White	10804-048	-	100
120 mL*	53 mm	Natural	10804-046	-	300

*Includes thiosulfate tablet for water testing.

