



Item Number B1194
Item Butylated Hydroxyanisole, Flakes, NF
CAS Number 25013-16-5
Molecular Formula $C_{11}H_{16}O_2$
Molecular Weight 180.25
Assay..... Min. 98.5%

MAXIMUM LIMITS

Residue on Ignition..... 0.01%
Heavy Metals.....0.001%
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	