



Certificate of Analysis

Teknova Controlled Document:CofA

DVC-Rev:	10674-02
Rev By:	JLA
Check By:	BT

Product Name / Specification:	DNA Suspension Buffer, pH 8.0. 250mL, RNase/DNase Tested, PCR Grade, Sterile.
-------------------------------	---

Version:
9.0

Part Number:	T0223
LotNumber:	T022312L2301
Date Manufactured:	11-Dec-2023
Date of Expiration:	11/2025

Date QC Released:	14-Dec-2023
-------------------	-------------

Batch Released By:

William Cummings, Senior
Manager, OA of OC

Appearance:	Clear, colorless liquid
Method of Sterilization:	0.2µm Sterile filtration
Storage Conditions:	Room temperature, 2-30°C.

QA Testing:	Result:	Expected:
Conductivity	571.2µS/cm	N/A
Density	0.99880 g/cm ³	N/A
DNase	Pass	Pass
PCR	Pass	Pass
pH	8.00	8.00 ± 0.10
RNase	Pass	Pass
Sterility	Pass	Pass