


AC-1272xx	CUSTOMER DRAWING	Page 1 of 1
Revision 03	Cryo.s™ 4 ml, Round Bottom, Starfoot Base	 greiner bio-one
	Item-No. 127 2xx	

NOT TO SCALE



Item No.	Description	A	B	C	D	E	F
127 261	Cryo.s™ 4 ml	Ø 9.8	Ø 12.4	Ø 13.1	76	83.3	14.5
127 263	Cryo.s™ 4 ml	Ø 9.8	Ø 12.4	Ø 13.1	76	83.3	14.5
127 277	Cryo.s™ 4 ml	Ø 9.8	Ø 12.4	Ø 13.1	76	83.3	14.5
127 278	Cryo.s™ 4 ml	Ø 9.8	Ø 12.4	Ø 13.1	76	83.3	14.5
127 279	Cryo.s™ 4 ml	Ø 9.8	Ø 12.4	Ø 13.1	76	83.3	14.5
127 280	Cryo.s™ 4 ml	Ø 9.8	Ø 12.4	Ø 13.1	76	83.3	14.5

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date September 03, 2009	Date September 03, 2009	Date September 04, 2009	
Date November 21, 2008	Name F. Falter	Name F. Sonntag	Name F. Sonntag	