

# adscope® Stethoscopes

## AFD Technology

Hear high- and low-frequency body sounds with Adjustable Frequency Diaphragms



### Lifetime Warranty and Free Parts Forever

*Electronics on the 658 are warranted for two years.*

#### CARDIOLOGY

#### CLINICIAN

#### PEDIATRIC



Model	658	600	601	602	606	615	608	603	612	619	618	614	604	605
Performance	10+	10	10	9.5	9.5	9.5	9.5	9	8	8	9	8.5	9	8.5
Main Feature	Electronic Amplification & Filter	AFD, Sculpted Chestpiece	AFD, Convertible	Traditional Cardiology	AFD, Ultra Lightweight	AFD, Sculpted Chestpiece	AFD, Convertible	AFD, Comfort	AFD, Lightweight	AFD, Ultra Lightweight	AFD, Patient Friendly	AFD, Sculpted Pediatric	AFD, Pediatric	Infant Size
Headset	Cardiology	Cardiology	Cardiology	Cardiology	Cardiology	Clinician	Clinician	Clinician	Clinician	Clinician	Clinician	Clinician	Clinician	Clinician
Chestpiece Material	Composite	Stainless	Stainless	Stainless	Aluminum	Zinc	Stainless	Stainless	Zinc	Aluminum	Composite	Zinc	Stainless	Stainless
Aural Tube Material	Stainless	Stainless <sup>1</sup>	Stainless <sup>1</sup>	Stainless	Aluminum	Stainless <sup>1</sup>	Stainless <sup>1</sup>	Stainless <sup>1</sup>	Stainless	Aluminum	Aluminum	Stainless	Stainless	Stainless
Length (in/cm)	30" / 76 cm	27" / 69 cm	28" / 71 cm	28" / 71 cm	28" / 71 cm	30" / 76 cm	30" / 76 cm	30" / 76 cm	30" / 76 cm	30" / 76 cm	30" / 76 cm	29.5" / 75 cm	29.5" / 75 cm	29.5" / 75 cm
Weight (oz/gr)	6 oz / 175 g	7.2 oz / 203 g	7.6 oz / 215 g	7.2 oz / 205 g	4.5 oz / 130 g	5.4 oz / 152 g	5.8 oz / 164 g	5.8 oz / 164 g	4.9 oz / 139 g	3.7 oz / 105 g	3 oz / 85 g	3.8 oz / 108 g	4.6 oz / 131 g	4.1 oz / 117 g
Chestpiece Finishes	N/A	Mirror, Gold Plated, Iridescent, Tactical	Brushed, Iridescent, Mirror, Rose Gold, Tactical	Brushed	Satin, Rose Gold, Tactical	Brushed, Tactical	Brushed, Copper, Iridescent, Mirror, Rose Gold, Tactical	Brushed, Gold-Plated, Copper, Iridescent, Rose Gold, Tactical	Brushed	Satin, Tactical, Champagne, Rose Gold, Smoke	N/A	Mirror	Brushed	Brushed
Extra Eartips	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	1 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus	2 pr. Adsoft Plus
Extra Diaphragm	No	No	No	Yes	No	No	No	No	No	No	No	No	No	Yes
Diaphragm Diameter	2"	1 3/4"	1 3/4" & 1 3/8"	1 3/4"	1 3/4"	1 3/4" Oval	1 3/4"	1 3/4"	1 3/4"	1 3/4"	1 1/2"	1"	1 3/8"	1 1/8"
Bell Diameter	N/A	N/A	1 1/2"	1 3/8"	1 3/8"	N/A	1 3/8"	1 3/8"	N/A	1 3/8"	N/A	N/A	1 1/8"	7/8"
Colors / Finishes	9	13	17	4	10	26	36	33	5	23	3	4	6	6
I.D. Tag Included	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Warranty	Electronics: Two Years Headset: Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime
COMPARABLE MODELS														
Littmann <sup>2</sup>	3100	Master Cardiology	Cardiology IV	N/A	N/A	Master Classic	Classic III	Classic II	Select	Lightweight	N/A	N/A	Classic II Pediatric	Classic II Infant
MDF <sup>2</sup>	N/A	N/A	797DD	797	N/A	N/A	777DT	777	N/A	747XP	N/A	N/A	777C	777I

<sup>1</sup> Some special finishes may have aluminum binaurals.

<sup>2</sup> Littmann is a registered trademark of the 3M Corporation. Master Cardiology, Cardiology IV, Master Classic, Classic II, and Classic III are trademarks of the 3M Corporation. MDF is a registered trademark of MDF Instruments.