



Quality Certificate

VWR® HDPE MULTIPURPOSE CONTAINERS

As per the manufacturer, the below product meets the following criteria:

| | |
|-----------------------------------|-------------------------------|
| North American Catalog No: | 89038-282 |
| Lot Number: | 43029 200802 |
| Description: | VWR Container Hdpe W/Lid Cs25 |
| Date of Manufacture: | August 02, 2020 |
| Country of Origin: | China |

Quality System Compliance

Products are manufactured under the ISO 9001:2015 & ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing:

Representative production samples are collected and inspected in accordance with current applicable standard protocols.

Dimensional Analysis: ASQ Z1.4-2003, General Insp. Level- II, AQL- 0.65. **Result:** Within tolerance

Product Specifications

This product is Class I category device as defined by the FDA in 21CFR Parts 862-892.

Material: **Container:** High-Density Polyethylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed
Lid: Blend of High-Density Polyethylene and Polyolefin Elastomer, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed

Expiration Date: For our non-sterile products, no expiration date (Shelf life) has been determined.

BSE/TSE Statement: All materials used and contact materials are of fully synthetic origin.

Non-Pyrogenic: N/A

BPA Statement: N/A

DEHP Statement: N/A

DNase & RNase Free: N/A

Cytotoxicity: N/A

Latex-free Statement: N/A

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 14, 2020