

## Certificate of Conformity

---

### Product Information

---

**Product Number:** 28903531

**Product Name:** illustra™ TempliPhi™ Sequence Resolver Kit (200 reactions)

**Lot Number:** 16976314

**Expiry Date:** 13 January 2021

The information contained in this certificate of conformity is GE Healthcare confidential and is subject to the use restrictions set forth in the Agreement between the parties.

### Contents

---

25640101A: Sample Buffer 1 x 900 µl Lot 16976314

25640101B: Reaction Buffer 1 x 900 µl Lot 16976314A

25640101C: Enzyme Mix 1 x 100 µl Lot 16976314

### Conformance & Quality systems statement

---

This is to certify that this product conforms to GE Healthcare Life Sciences specifications.

All products are manufactured via a Quality System certified to **BS EN ISO9001:2008** and tested in accordance with documented quality procedures and approved as a result of meeting the required specification

### Functional assay

---

A 20 µl amplification reaction of TempliPhi™ Sequence Resolver Kit must generate 0.1-0.5 µg of DNA in 4 hours and 1.4-2.4 µg in 18 hours from 10 ng of purified pUC19 DNA.

### Electronic signature

---

This document has been electronically produced and is valid without a signature.

Version AH