



CERTIFICATE OF ANALYSIS

Dehol[®]70

ISOPROPYL ALCOHOL SOLUTION

Lot # A06112002A

Manufactured: June 2020

Expiration: June 2022

Decon 8416 16 oz. Size

Meets Current USP Specifications

TEST	SPECIFICATION	UNIT	RESULT
Appearance	Clear, colorless liquid	N/A	Pass
Isopropyl alcohol, % v/v	68.0 – 72.0	%	70.6
Specific gravity at 20°C	0.872 - 0.883	N/A	0.877
Limit of non-volatile residues	< 5	mg	< 5
Acidity	<0.02	meq/50ml	0.01

The above material has been filtered through a Hydrophobic Cartridge filter, pore size 0.2 µm.

Reviewed and Approved by:

Date: June 2020

A handwritten signature in blue ink that reads 'Karen Wend'.

Karen Wend
Manager, Quality and Technical Services

Updated November 16, 2011- USP note

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.