

Quality Assurance Certificate

Microsart® Funnel 100

Order no. 16A07--10-----N

Use before 02 / 2025

Lot no. 211100487

This document is to verify that the designated product was manufactured by Sartorius Stedim according to a Quality Management System that is certified for compliance with current ISO 9001.

Sartorius Stedim hereby declares that the designated product was manufactured to meet the current specifications for these products. Sartorius Stedim also declares that the designated products have been manufactured from approved materials of construction and comply with design and manufacturing documentation.

Sartorius Stedim certifies that Microsart® Funnel are manufactured for use in the bacteriological analysis of potable waters, waste and natural waters and other liquids in accordance with procedures referenced in Standard Methods (recent edition) and laboratory standards for equipment and materials set forth by the U.S. Environmental Protection Agency.

- We herewith certify that the graduation marks on the Microsart® @filter 100ml for all manufactured batches have a tolerance at nominal fill capacity of $\pm 5\%$.
- Intermediate graduation marks have a tolerance of $\pm 10\%$ in accordance with ISO/DIN 7056 International Standards for Laboratory Plasticware Accuracy.


Specifications	Result
Visual inspection	Pass
Sealing test	Pass
Water flow test	Pass

Biological Properties
Sterility

Units have been sterilized with EO to a Safety Assurance Level (SAL) of 10^{-6}

14. April 2021

Date



Philip Stafford
Quality Assurance

Manufactured by Sartorius Stedim Lab Ltd.
GL 10 3UT Stonehouse, Gloucestershire, UK

Phone +44.1453.821972
Fax +44.1453.827928

www.sartorius.com