

Faster.
Safer.
More efficient harvesting.
Now it's all in the bag.



500 mL HarvestLine Bottle Liners

Safe, efficient workflows are as critical to your study as quality samples.

Now you can harvest high-quality samples in fewer steps from prep to cleanup.

With new HarvestLine sterile bottle liners, optimized speed, safety and pelleting efficiency are in the bag.

Optimize lab safety and efficiency

- Prevent contamination
- Reduce exposure to potentially hazardous biologics
- Reduce repetitive motions
- Prevent liquid spillage
- Reduce needed storage space

Accelerate workflows

- Achieve maximum sample yield
- Minimize non-value-added tasks e.g., bottle cleaning, autoclaving, scraping
- Expedite freeze/thaw lysing



It's all in the **bag**

500 mL Harvestline System Liners



Ordering Information

VWR Cat. No.	Description	Required Adapter
Liners		
76292-128	500 mL HarvestLine System Liner, qty 6 (non-sterile)	Material: Polyetherurethane, USP 23 Class VI
76292-130	500 mL HarvestLine System Liner, qty 6 (sterile w/peel pack)	
Bottles		
BK361691	Assembly, Polypropylene, 500 mL, 69 x 160 mm (qty. 6)	Required Adapter
BK361690	Assembly, Polycarbonate, 500 mL, 69 x 160 mm (qty. 6)	
Fixed Angle Rotors		
BK369687	JA-10 Rotor; 10,000 rpm; 17,700 x g	N/A
BK369681	JLA-10.500 - Biosafety	N/A
BK366754	JLA-9.1000 Rotor; 9,000 rpm; 16,800 x g	BK363683
BK363688	JLA-8.1000 Rotor; 8,000 rpm; 15,970 x g	BK363683
Swinging-Bucket Rotors		
BK339087	JS-5.2 Swinging-Bucket Rotor; 5,200 rpm; 6,840 x g	BK339109
76102-916	JS-4.750 Swinging-Bucket Rotor; 4,750 rpm; 5,250 x g	BK349945
BK362734	JS-4.3 Swinging-Bucket Rotor; 4,300 rpm; 4,220 x g	BK349945
BK339080	JS-4.2 Swinging-Bucket Rotor; 4200 rpm; 5,020 x g	BK339109
BK366695	JS-4.2A Swinging-Bucket Rotor; 4,200 rpm; 5,020 x g	BK339109
BK339086	JS-4.0 Swinging-Bucket Rotor; 4,000 rpm; 4,050 x g	BK339109
BK339081	JS-3.0 Swinging-Bucket Rotor; 3,000 rpm; 2,560 x g	BK339109
BK360581	GH-3.8 Swinging-Bucket Rotor; 3,750 rpm; 3,210 x g	BK349945
BK369702	SX-4750 Swinging-Bucket Rotor; 4,750 rpm; 5,250 x g	BK349945
BK369704	SX-4750A Swinging-Bucket Rotor; 4,750 rpm; 5,250 x g	BK349945
Accessories		
76292-136	Kit, Liner, 500 mL HarvestLine	
76294-188	– Adapters, HarvestLine 500 mL (qty 6)	
76294-190	– Base Spacers, HarvestLine 500/750 mL (qty 2)	
76294-198	– Filling Rack (stainless steel)	
76294-200	– Funnel (polyethylene)	

Visit vwr.com/beckman or contact your VWR Sales Representative for more information.

800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.

