

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Spinbar® Magnetic Stir Bars, Egg-Shaped
Country of Origin	Manufactured in United States
Date of Issue (yyyy-mm-dd)	2021-09-15

VWR Catalogue Number	Diameter	Length	For use with
58949-012	5 mm ($\frac{3}{16}$ ")	9.5 mm ($\frac{3}{8}$ ")	10 mL Flasks
58949-010	6 mm ($\frac{1}{4}$ ")	15 mm ($\frac{5}{8}$ ")	10 mL Flasks
58949-006	9.5 mm ($\frac{3}{8}$ ")	19 mm ($\frac{3}{4}$ ")	10, 25, 50, 100 mL Flasks
58949-196	12 mm ($\frac{1}{2}$ ")	25 mm (1")	100, 200 mL Flasks
58949-016	15.9 mm ($\frac{5}{8}$ ")	38.1 mm (1 $\frac{1}{2}$ ")	300, 500 mL Flasks
58949-200	15.9 mm ($\frac{5}{8}$ ")	32 mm (1 $\frac{1}{4}$ ")	300, 500 mL Flasks
58949-210	19 mm ($\frac{3}{4}$ ")	41 mm (1 $\frac{5}{8}$ ")	500, 1000 mL Flasks
58949-002	19 mm ($\frac{3}{4}$ ")	51 mm (2")	2000, 3000 mL Flasks
58949-243	19 mm ($\frac{3}{4}$ ")	64 mm (2 $\frac{1}{2}$ ")	3000, 5000 mL Flasks
58949-008	19 mm ($\frac{3}{4}$ ")	76 mm (3")	12000, 22000 mL Flasks

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

- Autoclavable @121° C
- Tolerance +/- .020"

Product Specifications

Material	Dupont White Teflon #807N X, Magnet, ALNICO 5 Teflon used to encapsulate the magnetic bar meets USP Class VI, compliance with FDA Regulations 21CFR177.1550, 21CFR177.1520. Above product does not contain or come in contact with lead, latex, melamine, Mercury, Formaldehyde or other contaminants.
Sterilization	Not applicable.
ATP Assay	Not applicable.
Melamine Statement	This product is free from melamine.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not applicable.
Storage Conditions	Store at ambient warehouse temperatures, low humidity, protected from direct sunlight.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.