

**Catalog Number:** BDH1129-4LP, BDH1129-4LG, BDH1129-19L, BDH1129-204L,  
CABDH1129-4LP, CABDH1129-4LG, CA BDH1129-19L, CABDH1129-204L,

**Catalog Description:** 60% n-Hexane

**Grade:** ACS


**Date Current Specification Set:** April 19 2006

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

### Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance	Clear
Water (K.F.):	Less than 0.010%
Residue:	Less than 0.001% (0.0010 g/145 mL)
GC Purity:	Greater than 98.5% n-hexane (n-C6) and C6 isomers
n-Hexane:	Greater than 60.0%
Color (APHA):	Less than 10
Water-soluble titrable acid:	Less than 0.3 µeq/g (Not more than 1.0 mL of 0.01N sodium hydroxide)
Sulfur compounds (as S):	Less than 0.005% (Absence of precipitate)
Thiophene:	Passes test (Acid layer should not be blue or green)

Signed on behalf of VWR:



Neal Fox  
Quality Control Supervisor