



Anachemia
A VWR Company

Product Specification

Material 48915-700
Material description Ionall, buffer solution

Additional information

Characteristics	Specifications
Appearance	Clear colourless liquid
pH (25°C)	5,0 --> 5,5
F (Fluoride)	Max. 0,10 ppm

Signature	Additional information
<p>We certify that this batch conforms to the specifications listed above.</p> <p>This document has been produced electronically and is valid without a signature.</p> <p>Michel Boisvert, Head of the Laboratory - Lachine VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3, QC, Canada</p>	<p>APHA 4500-F- C. for fluoride. Ion-Selective Electrode Method</p>