



# Certificate of Analysis

**Product:**
**Ficoll-Paque™ PLUS**
**Code Number:**

17144003

**Lot No:** 10298654

Test/Characteristic:	Limit:	Result:
<b>Sterility, F<sub>0</sub>-value</b> minutes	min. 15	Passed
<b>Density</b> g/ml, 20°C	1.076 - 1.078	1.077
<b>pH</b>	5.5 - 7.5	6.7
<b>UV / Visible spectra</b>	to pass test	Passed
<b>Sodium Diatrizoate</b> g/dl	8.1 - 9.9	9.1
<b>Colour; A = 405 nm</b> AU	max. 0.20	0.11
<b>Endotoxin activity</b> EU/ml	less than 0.12	<0.12

Note: Sterility is ensured by in-process control.

**Manufactured in compliance with our ISO 9001 certified quality management system.**

**No animal derived material is used in raw materials or in the manufacture of this product.**

**Approval Date (Year-Month-Day):** 2020-09-21

**Expiry Date (Year-Month):** 2023-09

**Manufacturing Date (Year-Month):** 2020-09

Cytiva and the Drop logo are trademarks of Global Life Sciences  
IP Holdco LLC or an affiliate

Cytiva Sweden AB  
Björkgatan 30  
SE-751 84 Uppsala, Sweden  
T + 46 (0)18 612 00 00  
Cytiva.com  
Reg.No. SE 55 61 08 1919 01

**Tests and limits according to AS 45-6006-81 Ed. AE**
**Issued (Year-Month-day): 2020-09-21 by Quality Assurance**
**Medin, Britta**

This document has been electronically produced and is valid without a signature.

DOC1103901 / 7