

## Quality Certificate

### VWR® SPINBAR® MAGNETIC STIR BARS, POLYGON WITHOUT PIVOT RING

As per the manufacturer, the below product meets the following criteria:

<b>North American Catalog No:</b>	58948-986
<b>Lot Number:</b>	J252441
<b>Description:</b>	VWR Bar Stir Magnetic Ptfе 2 1/2x1/2 In
<b>Date of Manufacture:</b>	November 14, 2019
<b>Product Shelf Life:</b>	Not Applicable
<b>Country of Origin:</b>	United States of America

### Quality Conformance

The product is in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation.

### Product Specifications

**Material:** Dupont White Teflon #807N-X, Magnet, ALNICO 5

**Dimensions:** 2 ½" x ½" (63.5mm x 12.7mm), +/- .020"

**Other Compliances:**

- Autoclavable at 121 °C
- Materials Meet FDA Regulations 21CFR177.1550, 21CFR177.1520
- Teflon used to completely encapsulate the magnetic bar meets US Pharmacopia Class VI
- Based on the information provided to us this product does not contain or come in contact with lead, latex or other contaminants.
- To the best of manufacturer's knowledge materials used do not contain or come in contact with any animal origin.
- All Items are carefully inspected to conform with ISO 9001:2008.
- There is no expiration for this product.

**Signed:**



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** December 30, 2019