

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Graduated Urine Centrifuge Tubes, Polystyrene
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-02-23

Tubes

VWR Catalogue Number	Description	Volume	Cap type	Max. rcf
21008-035	Bulk Packed Tubes	15 mL	Plug	1500 × g
21008-000	Tubes in Partitioned Boxes	15 mL	Plug	1500 × g

Caps

VWR Catalogue Number	Description
21008-022	LDPE Plug-Type Caps

Quality System Compliance

Products are manufactured under the [ISO 9001:2015 & ISO 13485:2016](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

This product is Class I category device as defined by the FDA in 21CFR Parts 862-892.

Material
Tube: PS (Polystyrene), Compliance with USP Class VI Color: Clear
Cap: LDPE (Low Density Polyethylene), Compliance with USP Class VI Colour: Natural

Sterilization Not applicable.

ATP Assay Not applicable.

DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993 standards for cytotoxicity and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is 0.05 EU/ml . (TAL Gel Clot Method).
Latex Statement	The product is Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Not applicable.
Storage Conditions	Store at room temperature.
Shelf-Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.