



## CERTIFICATE OF ANALYSIS

# CiDehol<sup>®</sup>70

### ISOPROPYL ALCOHOL SOLUTION

Lot # **A08122002D**

Manufactured August 2020

Expiration August 2022

Decon 8432 32 oz. Size

Meets Current USP Specifications

TEST	SPECIFICATION	UNIT	TEST RESULT
Appearance	Clear, colorless liquid	N/A	Pass
Isopropyl alcohol, % v/v	68.0 – 72.0	%	68.2
Specific gravity at 20°C	0.872 - 0.883	N/A	0.882
Limit of non-volatile residues	< 5	mg	< 5
Acidity	<0.02 meq/50ml	meq/50ml	<0.01

The above material has been filtered through a Hydrophobic Cartridge filter, pore size 0.2 µm.

Reviewed and Approved by:

Date: August 2020

A handwritten signature in blue ink, appearing to read 'Karen Wend'.

Karen Wend  
Manager, Quality and Technical Services

Updated November 16, 2011- USP note

**BSE/TSE Statement:** There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.

460 Glennie Circle • King of Prussia, PA 19406 • (610) 755-0800 • Fax (610) 270-8905