

New Products for Distribution
through VWR



Conical Competence

Eppendorf Conical Tubes 15 mL and 50 mL

Eppendorf Tubes offer a premium quality tube that meets the highest demands of any microbiology or molecular biology laboratory.

Eppendorf Conical Tubes 15 mL and 50 mL

- > Sterile, pyrogen-, DNase-, RNase-, and DNA-free
- > Each lot is checked and certified for quality and purity by an independent laboratory

Eppendorf Conical Tubes 15 mL and 50 mL DNA LoBind

- > Maximizes the recovery of nucleic acids by significantly reducing sample-surface binding
- > Ideal for sample preparation and long-term storage

Your local distributor:

Eppendorf Conical Tubes 15 mL and 50 mL




Innovative cap design

- > Pronounced inner lid ring ensures a tight seal and minimizes sample loss
- > Promotes user-friendly, ergonomic handling
- > Flattened edges on the lid support a stable, upright position—the lid won't roll away—minimizing the contamination risk
- > High contrast lid for easy label reading



High centrifugation stability

- > Up to 19,500 x g
- > Enables shorter centrifugation time to get your results faster
- > Ensures high safety for samples, user, and equipment



Premium material free of unnecessary additives

- > USP Class 6 raw material free of slip agents, plasticizers and biocides—substances which may leach into your sample causing deleterious downstream affects
- > High transparency for optimal sample visibility
- > Large labeling surface on the side for easy labeling of your solutions and samples

Ordering information

Description	VWR Catalog No.
Conical Tubes, 15 mL, sterile ¹ , bulk, case of 500 (10 bags x 50 tubes)	75875-238
Conical Tubes, 50 mL, sterile ¹ , bulk, case of 500 (20 bags x 25 tubes)	75874-294
Conical Tubes, 15 mL, PCR Clean ² DNA LoBind, case of 200 (4 bags x 50 tubes)	75874-296
Conical Tubes, 50 mL, PCR Clean ² DNA LoBind ² , case of 200 (8 bags x 25 tubes)	75874-298

¹ Certified sterile and free of pyrogens, DNases, RNases and DNA.

² Certified free of hDNA, DNases, RNases, and PCR inhibitors.

For ordering information

1.800.932.5000
vwr.com

Prices and product details are current when published, subject to change without notice.
 | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

0617 Lit. No. 050021W

For product and technical information

www.eppendorf.com · Email: info@eppendorf.com

In the U.S.: Eppendorf North America, Inc. 800-645-3050

ENA.C4.0207.A.US.17.VWR.TUB © 2017 Eppendorf AG. Eppendorf® is a registered trademark of Eppendorf AG. VWR, forms of VWR and the VWR logo and/or design are either registered trademarks, trademarks or service marks of VWR International, LLC in the United States and/or other countries. Product appearance, specifications, and/or prices are subject to change without notice.