

**Sterile Alcohol
Sterile Disinfectants
Sterile Water**

70% **Medical Device** **Research** **Sterile** **Process Equipment** **Alcohol** **70% EtOH** **Sterile Water**
Biotechnology **Wipes** **Disinfectants** **6% Hydrogen Peroxide** **Ready-to use** **Biotechnology**
Healthcare **IPA** **Laboratory** **Sprays** **Water** **Pharmaceutical** **70%** **Alcohol**
Pathogens **99% IPA** **EtOH** **Cleanroom** **Process Equipment** **IPA** **EtOH**
Water **Ready-to-use** **IPA** **Water 99% IPA**
Purified Water **6% Hydrogen Peroxide**



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com



Customer Service:
800-932-5000

STERILE ALCOHOL

CiDehol® ST

Sterile 70% Isopropyl Alcohol

CiDehol® ST is a ready-to-use 70% IPA solution in non-aerosol trigger spray and one gallon bottles. It is filtered to 0.2 µm, individually double-bagged and gamma-irradiated. Each lot is tested for bacterial endotoxin and sterility. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing. Made with WFI. EPA registered.



STERILE DISINFECTANTS & CLEANERS

Peroxigen™ ST

6% Sterile Hydrogen Peroxide

Sterile, double-bagged, gamma-irradiated 6% Hydrogen Peroxide solution made with Water-For-Injection (WFI) for use in cleanrooms and critical areas requiring a sterile cleaner. No residues, low-toxicity, non-flammable, and zero VOC profile help to assist in VOC compliance. Each trigger spray bottle is double-bagged with a sterile spray head and each case label contains a QR Code that is linked to a Lot Specific Document detailing parameters of QC, irradiation, and sterility testing. Not EPA registered.



Sanihol® ST

Sterile 70% Denatured Ethyl Alcohol

Sanihol® ST is a ready-to-use 70% Ethanol solution in non-aerosol trigger spray and one gallon bottles. It is filtered to 0.2 µm, individually double-bagged and gamma-irradiated. Each lot is tested for bacterial endotoxin and sterility. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing. Made with WFI. EPA registered.



CiDehol® ST 99

Sterile 99% Isopropyl Alcohol

CiDehol® ST 99 is a ready-to-use 99% IPA solution in non-aerosol trigger spray and one gallon bottles. It is filtered to 0.2 µm, individually double-bagged and gamma-irradiated. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing.



STERILE WATER




AquaPur™ ST

Purified Water USP Grade

AquaPur™ ST is sterile purified water for use in dilution of sterile disinfectants, manufacturing, general laboratory and cleanroom applications. Purified water is filtered to 0.2 µm, irradiated to a SAL 10⁻⁶ and tested for sterility and endotoxin using current USP method. Meets USP specs. Each case label contains a QR Code that is linked to a Lot Specific Document which details parameters of QC, irradiation, endotoxin and sterility testing. Not for injection, irrigation or infusion.



PRODUCT GUIDE

PRODUCT		SIZE	VWR CAT. #	CASE SIZE
CiDehol® ST 70% <i>Sterile 70% IPA</i>		12 oz. Trigger Spray	89176-254	16 x 12 oz.
		16 oz. Trigger Spray	89125-194	12 x 16 oz.
		32 oz. Trigger Spray	89125-196	12 x 32 oz.
		1 gal. Refill	89125-198	4 x 1 gal.
Sanihol® ST 70% <i>Sterile Denatured Ethanol</i>		16 oz. Trigger Spray	89125-218	12 x 16 oz.
		1 gal.	89125-220	4 x 1 gal.
CiDehol® ST 99% <i>Sterile IPA</i>		16 oz. Trigger Spray	89125-202	12 x 16 oz.
		1 gal.	89125-204	4 x 1 gal.
Peroxigen™ ST <i>Sterile 6% Hydrogen Peroxide</i>		16 oz. Trigger Spray	75853-810	12 x 16 oz.
AquaPur™ ST <i>Sterile Water</i>		1 gal.	71003-690	4 x 1 gal.



460 Glennie Circle
 King of Prussia, PA 19406
 (800) 332-6647
www.deconlabs.com



Customer Service:
 800-932-5000