

MATERIAL SAFETY DATA SHEET

HF scientific, inc.
3170 Metro Parkway
Fort Myers, FL 33916-7597
(239) 337-2116
Mon-Fri 8:00 am - 4:30 pm

24 Hour Emergency
Telephone Number
CHEMTREC
1-800-424-9300

SECTION 1 PRODUCT IDENTIFICATION

PRODUCT NAME: 4000 NTU STOCK SOLUTION - FORMAZIN KIT

PART NUMBER: 50040, 70914

CHEMICAL FAMILY: NA

SYNONYMS: NA

FORMULA: Hydrazine Sulfate <0.1%
Hexamethylenetetramine <10%
Formazin Polymer <1%
Water >88%

SECTION 2 HAZARDOUS INGREDIENTS

Other Limits

Hazardous Ingredient(s)	OSHA PEL	ACGIH TLV	
<u>Recommended</u>			
Hydrazine Sulfate mg/kg CAS# 10034-93-2	Not listed	Not Listed	LD ₅₀ 601 (oral rat)
Hexamethylenetetramine CAS# 100-97-0 rat)	Not Listed	Not Listed	LD ₅₀ 9200 mg/kg (intravenous
Formazin Polymer CAS# Not Applicable	Not Listed	Not Listed	Not Listed

SECTION 3 PHYSICAL DATA

Boiling Point(F): Approximately 212 ° F (100 °C)
Specific Gravity (H₂O = 1): Not found
Vapor Pressure (mm Hg): Not applicable
Melting Point: Not found
Vapor Density (Air = 1): Not found
Evaporation Rate (Butyl Acetate = 1): Not found
Solubility in Water: Fluid layer is miscible
Appearance and Odor: White suspension - little odor
Percent Volatile By Volume: Not found

SECTION 4 FIRE AND EXPLOSION HAZARD DATA

Flash Point: Not flammable

Flammable Limits:

LEL: Not Applicable **UEL:** Not Applicable

Extinguishing Media: Suitable to cause of fire

Special Fire Fighting Procedures: NIOSHA approved respirator and protective clothing.

Unusual Fire and Explosion Hazards: None found

SECTION 5 HEALTH EFFECTS AND FIRST AID

Route(s) of Entry: Inhalation? No Skin? Yes Ingestion? Yes

Health Hazards (Acute and Chronic): Acute and chronic for residual hydrazine sulfate

Carcinogenicity:

NTP? Yes, for hydrazine sulfate

IARC Monographs? Yes, for hydrazine sulfate

OSHA Regulated? Yes, for Hydrazine CAS # 302-01-02, 1.3 mg/m³ TWA

Hydrazine sulfate, which may exist in residual amounts in formazin is a

possible cancer hazard based on tests with laboratory animals.

Signs and Symptoms of Exposure: Irritation of the eyes, skin, and mucous membranes. Gastrointestinal pain if swallowed.

Medical Conditions Generally Aggravated by Exposure: Allergies or sensitivity to hexamethylenetetramine.

EMERGENCY AND FIRST AID PROCEDURES:

Eyes: Flush eyes for at least 15 minutes. See a physician.

Ingestion: Get medical help immediately.

Skin: Wash well with soap and water. If irritation develops, see a physician.

Inhalation: Remove to fresh air.

SECTION 6 REACTIVITY DATA

Stability: Stable under normal conditions of use and storage

Conditions to Avoid: None found

Incompatibility (Materials to Avoid): Not found

Hazardous Decomposition of Product: None

Hazardous Polymerization: Will not occur under normal conditions of use and storage

Conditions to Avoid: None found

SECTION 7 SAFE HANDLING AND USE

Steps to be Taken in Case Material is Released or Spilled: Absorb with earth, sand, or similar material. Flush area with water.

Waste Disposal Method: DISPOSE OF IN ACCORDANCE THE APPLICABLE FEDERAL, STATE AND LOCAL ENVIRONMENTAL REGULATIONS.

Precautions to be Taken in Handling And Storage: Normal precautions satisfactory. Avoid skin and eye contact.

Other Precautions: Do not pipette by mouth

SECTION 8 CONTROL MEASURES

Respiratory Protection: None required in normal use

Ventilation: General ventilation is satisfactory for normal use

Skin Protection: Rubber gloves

Eye Protection: Splash goggles

Other Protective Clothing or Equipment: Lab coat suggested. Safety shower and eye wash.

Work/Hygienic Practices: Wash well after handling.

NO GUARANTEE IS MADE AS TO THE ACCURACY OF ANY DATA OR STATEMENT CONTAINED HEREIN. WHILE THIS MATERIAL IS FURNISHED IN GOOD FAITH, NO WARRANTY EXPRESSED OR IMPLIED, OR MERCHANTABILITY, FITNESS OR OTHERWISE IS MADE. THIS MATERIAL IS OFFERED ONLY FOR YOUR CONSIDERATION, INVESTIGATION AND VERIFICATION AND HF SCIENTIFIC, INC., INCLUDING ITS DIVISIONS, AFFILIATES AND SUBSIDIARIES SHALL NOT IN ANY EVENT BE LIABLE FOR SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH ITS PUBLICATION. LIKEWISE, NO STATEMENT MADE HEREIN SHALL BE CONSTRUED AS A PERMISSION OR RECOMMENDATION FOR USE OF ANY PRODUCTION IN A MANNER THAT MIGHT INFRINGE EXISTING PATENT.